Diabetic Ulcers

THE V.A.C.®
VACUUM ASSISTED CLOSURE™

Clinical Summary Sheet

Comparison of Negative Pressure Wound Therapy using V.A.C.® with advanced moist wound therapy in the treatment of diabetic foot ulcers. A Multicenter Randomised Controlled Trial.

Authors: Blume P.A., Walters J., Payne W., Ayala J., Lantis J.

Journal: *Diabetes Care*, Vol 31; No 4; 631-636; 2008

Background

- A disabling complication of diabetes is foot ulcer development which leads to nonhealing chronic wounds that are difficult to treat
- Diabetic foot ulcers (DFU) are a significant risk factor for nontraumatic foot amputations in individuals with diabetes.

Purpose

To evaluate the safety and efficacy of V.A.C.® Therapy compared with advanced moist wound care (AMWC) for the treatment of DFUs.

Methods

- Concealment of randomisation
- Off loading therapy
- Independent wound assessments

Main entry criteria

≥ 18 year old

Ulcer Wagner stage 2 or 3

Calcaneal, dorsal or plantar foot ulcer ≥ 2cm²

Adequate blood circulation (perfusion)

No active Charcot disease


No ulcer malignancy

No untreated osteomyelitis

No uncontrolled hyperglycemia

No normothermic or hyperbaric oxygen therapy within 30 days of study start


No use of corticosteroids, immunosuppressive medication or chemotherapy within 30 days of study start


Clinical Summary Sheet

Results and Conclusions


V.A.C.[®] Therapy leads to greater wound area reduction


V.A.C.[®] Therapy leads to higher rate of wound closure


V.A.C.[®] Therapy reduces time to wound healing

V.A.C.[®] Therapy shows faster wound bed preparation

V.A.C.[®] Therapy reduces the rate of secondary amputations

*The duration of the control treatment could not be estimated (Kaplan-Meier).

Clinical and Economical Summary Sheet

Resource utilisation and economic costs of care based on a randomised trial of vacuum-assisted closure therapy in the treatment of diabetic foot wounds.

Authors: Apelqvist J., Armstrong D.G., Lavery L.A., Boulton A.J.M
Journal: *The American Journal of Surgery*, 2008 Jun;195(6):782-8

Background

- Direct treatment costs for diabetic foot complications are high, exceeding \$30,000. The costs for performing a single amputation are approximately \$34,000.
- These costs do not account for increased utilisation of inpatient services, prosthetic/rehabilitation, or home care/social services.
- Societal costs are especially high in cases of amputations and ulcers that fail to heal as patients require prolonged hospitalisation, rehabilitation, need home nursing care, lose work production (for the disabled patient and his or her family).
- Important clinical factors related to short-term costs include healing rate and the presence and severity of foot infections and ischemia.

Purpose

To estimate and compare resource utilisation and direct economic costs of treatment with V.A.C.[®] Therapy in patients with severe postoperative diabetic foot wounds.

Methods

- Retrospective analysis of patients enrolled in a randomised clinical trial (Armstrong, Lancet 2005), using payer and societal perspectives.
- RCT with 162 patients with partial diabetic foot amputation of 16-week duration comparing V.A.C.[®] Therapy and Advanced Moist Wound Care (AMWC).
- Patients fulfilling a minimum of 8 weeks participation in the study period were included when comparing and analysing resource utilisation and economic calculations.


Resource Utilisation and Cost Calculation	
Inpatient care	<ul style="list-style-type: none">• Average cost per treatment day including resource use• Costs of surgical procedures (flaps, grafts, amputations, surgical debridements)
Antimicrobial treatment	<ul style="list-style-type: none">• Costs for oral and systemic antibiotics
Outpatient treatment visits	<ul style="list-style-type: none">• Number of outpatient clinic visits• Number of visits to patient's home• Estimated cost for caregiver time• Estimated cost of clinic visit
Topical wound treatment	<ul style="list-style-type: none">• Actual number of dressing changes• Estimated material costs (primary dressing material)• Estimated time consumption for dressing changes• Staff who performed dressing changes (cost/hour)

Clinical and Economical Summary Sheet


Results and Conclusions

- A higher proportion of wounds treated with V.A.C.® Therapy healed compared to AMWC (55.8% vs 38.8%, $p=0.040$).
- 89.1% of patients were treated in homecare setting after acute care intervention.
- V.A.C.® Therapy is a cost-effective intervention for patients with complex postoperative wounds.
- Treatment with V.A.C.® Therapy resulted in a greater proportion of patients obtaining wound healing at a lower overall cost of care when compared to AMWC.


V.A.C.® Therapy reduces the number of procedures


V.A.C.® Therapy reduces the number of outpatient visits


V.A.C.® Therapy reduces the number of dressing changes

Lower average direct treatment costs for V.A.C.® Therapy

Lower average total cost of treatment for V.A.C.® Therapy

VAC-THERAPY AND LASER-INDUCED FLUORESCENCE OF INDOCYANINE-GREEN (IC-VIEW), AN ASSESMENT OF WOUND PERFUSION IN DIABETIC FOOT SYNDROME

Zöch, Gerald MD,

Surgical Department, SMZO–Donauspital-Vienna Langobardenstrasse 122. A 1220 Vienna, Austria

Aim of the study was the use of IC-View for monitoring the effect of VAC-therapy on perfusion in diabetic wounds and the surrounding skin. We present data of ten patients with diabetic foot syndrome.

Methods:

IC-VIEW is a near-infrared video system, which visualises and quantifies tissue perfusion. The measurement is based on the fluorescence of the near-infrared absorbing and fluorescent agent indocyanine green (ICG). A measurement consists of following actions: Positioning the IC-VIEW in a distance of 50 cm from the tissue of interest, iv-injection of 0,5 mg/kg ICG, recording of ICG fluorescence for 1-8 minutes, quantification of fluorescence signal using the IC-CALC software and automated calculation of perfusion (%). Data for the improvement (%) of perfusion for the wound and the surrounding skin are shown.

Seven men and three women with age from 56 to 86 (mean 68) years were treated with VAC therapy for 7 to 24 (mean 12) days. The perfusion index was calculated for the wound itself and for the surrounding skin before and at the end of VAC-therapy.



Results:

The perfusion of the wound showed an increase under VAC-therapy of 31% (6-62,5%). The perfusion index calculated for the skin around the defect achieved an increase of 14% (1,5-38%). All defects were covered with meshed skin grafts and healed.

Conclusion:

IC-View enables to quantify the improvement of perfusion and granulation tissue in diabetic wounds due to VAC therapy.

Age/sex	VAC-therapy (days)	perfusion (%) of the skin			perfusion (%) of the wound		
		before	after	improvement	before	after	improvement
01 B.E./ 86 m	9	55	71	22	21	24	12,5
02 S.J./ 65 m	9	39	44	10	30	56	46
03 K.F./ 73 m	7	44	51	14	31	40	22,5
04 T.I./ 63 f	24	44	52	15	32	34	6
05 T.S./ 58 f	14	41	55	25	52	00	48
06 S. E./ 56 f	20	42	80	38	47	89	47
07 G.J./ 60 m	9	72	73	1,5	21	56	62,5
08 H.E./ 72 m	7	62	65	5	54	76	29
09 Z. F./ 76 m	13	40	42	6	48	51	6
10 T.J./ 71 m	7	48	51	6	35	42	31
mean 68	12	49,7	58,4	14,0	37,1	56,8	31

Topical Negative Pressure: Management of Complex Diabetic Foot Wounds

DG ARMSTRONG, AJM BOULTON & PE BANWELL

SERIES EDITOR: **PE Banwell FRCS(PLAST)**

CONTENTS

Introduction	2
Epidemiology	3
Etiology of Foot Ulceration in Diabetes	4
Principles of Management of Diabetic Foot Ulcers	5
Current Dressing Regimes for Diabetic Foot Ulcers	6
Topical Negative Pressure Therapy Theory	7
Topical Negative Pressure Components	8
TNP Dressing Application	9
Topical Negative Pressure Therapy: Mechanism of Action	11
Use of TNP in Chronic Wounds	14
Specific Considerations for the Use of TNP in Wounds of the Diabetic Foot	15
Conclusion	18
Bibliography	19

Topical Negative Pressure: Management of Complex Diabetic Foot Wounds

DG ARMSTRONG, AJM BOULTON & PE BANWELL

SERIES EDITOR: **PE Banwell** FRCS(PLAST)

CONTENTS

Introduction	2
Epidemiology	3
Etiology of Foot Ulceration in Diabetes	4
Principles of Management of Diabetic Foot Ulcers	5
Current Dressing Regimes for Diabetic Foot Ulcers	6
Topical Negative Pressure Therapy Theory	7
Topical Negative Pressure Components	8
TNP Dressing Application	9
Topical Negative Pressure Therapy: Mechanism of Action	11
Use of TNP in Chronic Wounds	14
Specific Considerations for the Use of TNP in Wounds of the Diabetic Foot	15
Conclusion	18
Bibliography	19

Introduction

Foot ulceration represents an all too frequent complication of diabetes, with neuropathy and peripheral vascular disease being important contributory factors in the pathogenesis of such ulceration.¹ However, it is not an inevitable consequence of having diabetes that a certain proportion of patients will develop foot lesions. Elliot Joslin, the famous Boston diabetes physician, realized this almost 70 years ago when, after observing many clinical cases of diabetic foot disease, he remarked "diabetic gangrene is not heaven sent, but earth born".² Foot ulceration therefore invariably occurs in certain individuals as a consequence of an interaction between specific pathologies of the lower limb and environmental hazards.

One of the most important contributory factors to the development of foot ulceration is the loss of pain, pressure, and joint position sensation that occurs as a consequence of diabetic neuropathy. The traditional medical model of disease suggests that a patient goes to the doctor because of symptoms, treatment is prescribed and then the patient recovers. Because this does not apply in the case of insensitive feet, it is often difficult for doctors to comprehend the diabetic foot syndrome: they find it difficult to take the initiative and look for early lesions or warning signs of imminent break-down. Many doctors regard such patients as being stupid — how can a sensible person walk on a swollen red foot with an active foot ulcer? What we need to realize is that an insensitive foot is not only painless, but also does not feel as if it belongs to the individual.³ Our task is to identify patients at risk of ulceration and help them to cope with this health state by avoiding exposure to environmental hazards that may result in ulceration.

Figure 1: A diabetic foot ulcer

Epidemiology

As stated above, foot ulceration is much more common in diabetic patients with neuropathy and vascular disease: the annual incidence rises from less than 1% in those without neuropathy for example, to more than 7% in those with established neuropathy.^{4,5} Globally, the diabetic foot remains a major medical, social, and economic problem that is seen in every continent.⁶ In developing countries, peripheral neuropathy is the most important risk factor for foot ulceration, and ulcers are frequently associated with quite severe limb sepsis in most of these countries.⁷

Naturally there are some regional differences between countries, for example the prevalence of risk factors for neuropathy appears to be slightly lower in Australia than in the United Kingdom.⁸ In the United States, it is recognized that there are differences in the epidemiology of diabetic foot problems according to ethnic background. Native Americans and Hispanics, for example, appeared to have a greater risk of lower extremity amputations amongst a population of diabetic veterans.⁹

In summary, diabetic foot ulceration and risk of lower extremity amputation remain all too common, but the study of the epidemiology of these problems has been clouded by a lack of agreement of a diagnostic criteria, and by variation in subject selection methods. In most populations, it must be assumed that at least 50% of older type 2 persons with diabetes have identifiable risk factors for foot ulceration.


Figure 2: General incidence of diabetes complications (adapted from Armstrong DG, Lipsky BA, 'Advances in the treatment of diabetic foot infections', *Diabetes Technologies & Therapeutics*, 2004, in press)

Current Dressing Regimes for Diabetic Foot Ulcers

There are a multitude of dressings and devices designed to aid in healing of the diabetic foot. While many have proven efficacious, others have been less so. Much of the reasons associated with less than desirable results relates not to the dressing or device itself, but rather to the attention given to the wound. Diabetic foot wounds must be effectively debrided of non-viable tissue and receive effective pressure off-loading in order to heal. Any modality, including TNP, will not live up to its full therapeutic potential without appropriate attention to debridement and pressure relief.¹³⁻¹⁵ With this caveat in place, there are a bewildering number of frequently employed dressings and techniques in diabetic foot care. They fall into numerous categories, which are outlined in Table I. The central tenet to all of the modalities listed below fall into the categories of infection control, drainage control, and control of the wound's microenvironment. Some modalities, such as topical negative pressure, are able to potentially influence all of these categories of control.

- | | |
|--|---|
| 1. Antimicrobial/Silver-containing Dressings | 15. Growth Factors |
| 2. Alginate Dressings | 16. Bioengineered Tissue Products |
| 3. Biosynthetic/Hyaluronan-based Dressings | 17. Hydrocolloid Dressings |
| 4. Cleansers | 18. Hydrogel Dressings (see also Impregnated Gauze Dressings) |
| 5. Collagen Dressings | 19. Transparent Film Dressings |
| 6. Composite Dressings | 20. Wound Fillers: Pastes, Powders, Beads, Strands |
| 7. Compression Bandages/Wraps | 21. Adhesive Skin Closures |
| 8. Conforming/Wrapping Bandages | 22. Topical Antibiotic Preparations |
| 9. Contact Layers | 23. Antimicrobials |
| 10. Hydrating Creams/Ointments/Hydrogels | 24. Moisture Barrier Ointments/Creams/Skin Protectant Pastes |
| 11. Devices | 25. Moisturizers |
| 12. Enzymes/ Debriding Agents | 26. Ointments |
| 13. Foam Dressings | 27. Biosurgical Modalities (Larvae) |
| 14. Gauze Dressings (see also Composite Dressings) | |

Table I: Partial list of categories of dressings and advanced wound healing modalities used in the diabetic foot¹⁶

Topical Negative Pressure Therapy Theory

Topical negative pressure therapy (TNP)¹⁷ has emerged as a powerful, non-pharmacological tool that can manipulate the wound healing process.¹⁸ This technique, based upon simple principles but delivering high-tech therapy, involves the application of a vacuum force across the wound surface, using a foam dressing interface. The physical forces generated at the foam-wound interface (via mechanical transduction) lead to a variety of changes in the wound environment, positively enhancing the healing process.

Numerous studies have now reported the successful use of topical negative pressure, known commercially as Vacuum-Assisted Closure™ (KCI, San Antonio, TX), in a clinical setting.^{19–26} Primary indications for use have recently been recommended¹⁸ and include the following: upper and lower limb trauma including thermal injury, infected wounds, sternal dehiscence, abdominal dehiscence, skin graft fixation, wound bed preparation, pressure ulcers, venous leg ulcers and diabetic ulcers. Furthermore, recent experience has shown useful applications in oral and maxillofacial surgery, gynaecological problems, necrotizing fasciitis, extravasation injuries and reducing donor site morbidity.¹⁸ However, on-going clinical trials must confirm the clinical and economic efficacy of this treatment.


Figure 5: TNP dressing *in situ*

- Traumatic wounds – upper/lower limb, burns
- Infected wounds
- Sternal dehiscence
- Abdominal dehiscence
- Skin graft fixation
- Wound bed preparation & exudate management
- Pressure ulcers – sacral, trochanteric, ischial
- Leg ulcers – venous, diabetic

Table II: Primary indications for using topical negative pressure therapy in clinical practice

Topical Negative Pressure Components

Several components are required for optimum delivery of topical negative pressure therapy: a foam interface dressing, a hermetically sealed drape, evacuation tubing, a wound effluent collection device and a microprocessor-controlled vacuum pump.

Specially-designed open pore, reticulated foam dressing should be used to cover the wound. These may be made of either polyurethane ether (PU) – 'black' foam, or polyvinyl alcohol (PVA) – 'white' foam. Each has individual characteristics as listed in Figure 10. Either type may be used effectively although the authors predominantly use PU foam in their practice.

An adhesive drape must be used to cover the foam dressing in order to create a closed, controlled wound environment. Often it may be easier to cut this drape into strips to facilitate securing the foam in place.

The evacuation tube is attached to the closed dressing by means of a so-called TRAC™ pad (KCI, Kidlington, UK) (TRAC = Therapeutic Regulated Accurate Care). This is a highly-sophisticated technology which allows a constant feedback from the wound site to the microprocessor-controlled vacuum pump V.A.C.® ATS™ (KCI, Kidlington, UK) and ensures that a constant negative pressure is delivered to the wound. Target pressures and delivery regimes are controlled by the vacuum pump which can be set to a variety of pressures using either continuous or intermittent cycles.

Figure 6: Microprocesor-controlled V.A.C.® ATS™ therapy unit (vacuum pump)


Figure 7: Foam dressings



Figure 8: TRAC™ technology


Figure 9: Dressing *in situ*

Characteristics

- Material: polyurethane foam
- Reticulated – coarse, open pore structure
- Hydrophobic foam
- Compression load deflection (firmness) 3kPa
- Tensile strength 160 kPa

Characteristics

- Material: polyvinyl alcohol foam, pre-moistened with sterile water
- Non-reticulated – denser, less open pore structure
- Hydrophilic foam
- Compression load deflection (firmness) 7 kPa
- Tensile strength 275 kPa

Figure 10: Differential characteristics of (left) polyurethane (PU) foam and (right) polyvinyl alcohol (PVA) dressing

TNP Dressing Application

Optimum results with topical negative pressure therapy can only be achieved if there is meticulous attention to detail regarding the dressing application. This is especially of paramount importance in dressing complex diabetic wounds.

The wound should first be debrided of all non-viable tissue and subsequently prepped using standard Institutional wound care regimes. This may occur in the operating room, on the ward, or in the community and should be considered one of the most important parts of the dressing regime.

Attention should also focus on peri-wound skin. It should be cleaned and dried. Barrier protectants, eg. Cavilon (3M) may be used as well as tincture of benzoin or friars balsam to ensure adhesiveness of the transparent drape. Additionally, many clinicians will also use hydrocolloid flanges to protect surround skin or augment the seal. This may also be useful to protect smaller wounds (<4 cm) from direct pressure effects of the TRAC™ pad.²⁷

Foam dressings (either PU or PVA) should be then be cut to size and placed within the wound. Purpose-made, indication-specific and anatomical dressings are now also being designed although certain techniques have been employed using existing foams.²⁸ There has been a trend to use an interposed dressing between the foam and the wound but this may dissipate the mechanical transduction effects at a given pressure and minimize granulation tissue formation. One of the rationales for using these dressings is to reduce the discomfort of removing the foam at dressing changes. However, in neuropathic conditions it would seem sensible to defer using them.

The foam should then be sealed with the adherent drape making sure additional shear stresses are not placed on surrounding skin causing traction injuries. Once sealed, a small hole is made in the drape and the TRAC™ pad applied before connecting up to the microprocessor-controlled vacuum pump. Every effort should then be made that a good seal has been applied and that the connectors are correctly in place. Adaptations of this technique may be employed which are especially useful for plantar wounds (see pages 16 and 17).

A standard protocol is currently employed for diabetic wound regimes: continuous pressures of 125 mmHg are routinely used although this may be increased or decreased according to wound response and patient comfort. It is recommended that dressing changes are performed every 48 hours or sooner in infected wounds. However, with increasing clinical experience and in non-infected or minimal-exudative wounds, dressing changes may be increased to fit in with institutional dressing regimes.

Certain wounds may also be suitable for instillation techniques using topical antiseptics and antibiotics.²⁹

Figure 11: A leg ulcer in a diabetic patient treated with TNP. Note the protective hydrocolloid flange

Topical Negative Pressure Therapy: Mechanism of Action

At a macroscopic level TNP can be regarded as a method, which combines the benefits of both closed and open wound management. Closed treatment by approximation of skin edges protects the wound from bacterial invasion and maintains moisture, pH, gaseous concentrations and temperature at a more optimal level for healing. Open treatment, however, avoids tension, allows for the release of toxic products thereby reducing the risk of systemic sepsis in the presence of contamination and enables regular monitoring and assessment of a complex wound. The unique features of TNP therapy, which are thought to contribute to an optimised wound environment include:

1. Removal of excess interstitial fluid:

There exist two mechanisms by which removal of excess interstitial fluid is thought to accelerate healing. The first is based on the observation that fluid from chronic wounds is known to contain cytokines, collagenases and elastases whose accumulation inhibit the proliferation of fibroblasts, keratinocytes and endothelial cells.³⁰⁻³³ It therefore follows that removal of such substances might serve to accelerate the wound repair process. Secondly, the localised oedema present in both acute and chronic wounds causes an increase in tissue turgor, which compresses surrounding microvascular and lymphatic channels. Diminution of interstitial fluid by TNP is thought to decompress small blood vessels thereby restoring their flow.^{34,35}

2. Delivery of mechanical stress to the wound:

This mechanism has been compared to those at play in the techniques of Ilizarov bone distraction³⁶ and tissue expansion and has led to TNP being likened to "reverse tissue expansion".¹⁷ At a cellular level, the mechanisms believed to be involved stem from Thoma's postulation that the application of mechanical stress results in angiogenesis and tissue growth.³⁷ This phenomenon appears to be controlled through cytoskeletal changes, which stimulate the release of intracellular second messengers capable of stimulating cell growth.³⁸⁻⁴⁰

3. Provision of a moist wound environment:

By virtue of its sealed system, TNP therapy creates a moist wound environment, which is thought to contribute further to accelerated healing. This concept, which came to light following Winter's classic experiments in the 1960s, was based on his observations that moist wounds were seen to re-epithelialise significantly more rapidly than those left exposed to air.^{41,42} His findings have been supported by numerous other investigators and have led to the current acceptance that moist wound healing is advantageous to the wound repair process.⁴³⁻⁴⁵

A full understanding of the exact microscopic mechanisms of action by which TNP exerts its beneficial effects are currently unknown but substantive evidence exists to support the following postulated mechanisms:

(a) Stimulation of angiogenesis and local blood flow:

This is realised by decompressing small blood vessels and stimulating the proliferation and differentiation of endothelial cells. Using Doppler probes in an animal model, Morykwas and Argenta demonstrated a four-fold increase in blood flow when 125 mmHg of negative pressure was applied to full-thickness porcine wounds.²⁰ This increase was maintained for 5–7 minutes when using continuous delivery and could be reproduced by intermittent periods of therapy. For the purposes of maximal blood flow the optimal treatment protocol was described as being intermittent delivery of 125 mmHg negative pressure for 5-minute periods with 2-minute intervals. Josephs *et al.* demonstrated enhanced angiogenesis histologically from biopsies of wounds treated with TNP⁴⁶ whilst Walgenbach and co-workers were able to demonstrate increased capillary formation using histochemical techniques to detect endothelial cell markers CD31 and CD105.⁴⁷

(b) Reduction in bacterial count:

Increased blood flow serves to increase tissue oxygen concentrations thereby inhibiting the growth of anaerobic organisms and providing nutrition for the oxydative bursts of cells of the immune system. Gouttefangeas and colleagues have demonstrated infiltration of the foam by granulocytes and mononuclear cells.⁴⁸ Neutrophils are known to form the first line of defense against wound infections,⁴⁹ however, in contrast to macrophages and T lymphocytes they are not essential for wound healing.^{50–52} Although there is no evidence to suggest that a higher concentration of such cells exists in a wound treated with TNP, the identification of their presence suggests a possible mechanism by which bacterial count is reduced as demonstrated by Morykwas and Argenta following challenge by inoculation with human isolates of *Staphylococcus aureus* or *Staphylococcus epidermidis*.²⁰


Figure 12: Mechanism of action (diagrammatic)

(c) Mitogenesis of cells involved in the wound repair process:

This proliferative response is thought to be due to a combination of mechanical force, moist environment and removal of inhibitory factors and is manifested by an increase in volume of granulation tissue which has been observed in both acute²⁰ and ischaemic animal models⁵³ of wound healing by as much as 63% with continuous TNP and over 100% with intermittent TNP. Mullner *et al.* conclude that TNP is capable of stimulating the production of granulation tissue to such an extent that musculocutaneous flap coverage was averted.²³


Figure 13: Foam/wound interaction illustrating granulation tissue formation

- **Removal of excess interstitial fluid**
- **Delivery of mechanical stress to the wound**
- **Provision of a moist wound environment:**
 - Stimulation of angiogenesis and local blood flow
 - Reduction in bacterial count
 - Mitogenesis of cells involved in the wound repair process

Table III: Mechanisms of action of TNP

(c) Mitogenesis of cells involved in the wound repair process:

This proliferative response is thought to be due to a combination of mechanical force, moist environment and removal of inhibitory factors and is manifested by an increase in volume of granulation tissue which has been observed in both acute²⁰ and ischaemic animal models⁵³ of wound healing by as much as 63% with continuous TNP and over 100% with intermittent TNP. Mullner *et al.* conclude that TNP is capable of stimulating the production of granulation tissue to such an extent that musculocutaneous flap coverage was averted.²³


Figure 13: Foam/wound interaction illustrating granulation tissue formation

- **Removal of excess interstitial fluid**
- **Delivery of mechanical stress to the wound**
- **Provision of a moist wound environment:**
 - Stimulation of angiogenesis and local blood flow
 - Reduction in bacterial count
 - Mitogenesis of cells involved in the wound repair process

Table III: Mechanisms of action of TNP

Use of TNP in Chronic Wounds

When using TNP in the management of chronic wounds a protocol of 125 mmHg continuous pressure has been recommended. Using this delivery strategy Joseph and colleagues⁴⁶ assessed 24 patients with 36 chronic wounds of mixed aetiology, which had been present for at least one month. Half of the wounds were treated with TNP and the other half with wet to moist saline dressings, which were changed three times a day. The change in wound volume when assessed at 6 weeks was 78% in the TNP-treated group and 30% in controls with a reduced rate of complications in the TNP-group (17% compared with 44%). Philbeck *et al.*⁵⁴ attempted to assess the cost implication of using TNP in the management of chronic wounds by comparing their cohort of individuals treated with low air loss beds and TNP with Ferrell's study⁵⁵ in which similar wounds were treated with low air loss beds and moist dressings. Although the average cost of hiring a V.A.C.[®] is \$100/day and the total cost of treating wounds with TNP per day was \$150 compared with \$95 overall costs were 38% less due to the more rapid (61% faster) rates of healing in this group.

In the literature review reported by Evans *et al.*⁵⁶ the authors describe a study in which 10 patients with diabetic foot ulceration of at least one month's duration were randomised to receive either TNP (125 mmHg continuously for 48 hours followed by intermittent therapy) or saline moistened gauze changed twice daily until definitive closure was achieved whether by secondary intention, skin grafting or flap coverage. Ulcers treated with TNP had an average period of 23 days until complete closure and a decrease in surface area of 28% compared with 43 days and 9% in the control group.^{57,58} Scheufler and colleagues⁵⁹ describe a successful two-phase TNP treatment protocol used on 19 lower leg wounds each of diabetic, ischaemic and traumatic aetiology. Mean surface area prior to débridement was 24 cm² and all defects were full-thickness with exposed tendon or deep muscle fascia in the base of the wound, which had been unresponsive to prior conventional therapy. Their protocol consisted of surgical débridement with intravenous antibiotics as indicated and intermittent TNP at 125 mmHg (3 minutes on, 2 minutes off) for between 1 and 4 weeks until a complete layer of granulation tissue was observed. In the second phase of treatment, gentle curettage of the granulation tissue was followed by skin grafting and TNP therapy at 125 mmHg continuous pressure for a further 4–5 days. The authors achieved a 90–100% take in 17 out of 19 patients with closure achieved within 3–5 weeks. Treatment failure in the remaining 2 patients was due to infection and both responded to a repeat of the protocol.

Figure 14: Non-healing leg ulcer


Figure 15: TNP stimulated granulation tissue prior to grafting. 100% graft take

Specific Considerations for the Use of TNP in Wounds of the Diabetic Foot

As mentioned in the preceding paragraphs, the central tenet to treatment of the non-infected, non-ischemic diabetic foot wound is appropriate debridement and offloading. Not adhering to these principles may be an explanatory variable for failure of therapy.^{13,14} The patient with a diabetic foot wound is not “protected” by pain, and therefore may not perceive many warning signs of skin breakdown.^{3,60} True appreciation of and accommodation for this neuropathic state is an essential component to superior wound management.

Many patients, however, have very deep, complex plantar diabetic foot wounds that are not initially amenable to simple offloading techniques. If there is not an appropriate plan for dead space management, these deep wounds are at higher risk for infection and amputation.⁶¹ Previous treatment for these large wounds has involved extended periods of complete non-weightbearing, frequent dressing changes, and convalescence in hospital or long-term care facilities. Topical negative pressure has shown some significant utility in this patient population.^{20,58,62–64} However, weight bearing on a TNP dressing has, in the past, been quite difficult because the tubing (contiguous with the polyurethane sponge) generally exits plantarly. Any weight placed on this construct could cause tissue damage.

In other parts of the body, clinicians have applied the TNP dressing using a technique known as “bridging” (Figures 16 to 19)^{28,65} as a strategy to connect anatomically distant wounds to a single vacuum unit. We have modified this technique to allow for weight bearing while still using TNP. This involves applying sponges from an area of potentially high pressure and/or repetitive stress to an area that would allow for a more simple and safe connection. To be specific, the sponge is generally run from the site of the plantar diabetic foot wound to the dorsum of the foot where the TNP tubing is connected to a portable vacuum unit. The patient is then placed into an offloading device such as a removable cast walker, which can be converted into an “instant total contact cast” (Figures 16–19).^{13,14,66} We have found that this allows a patient to perform many unaided ambulatory activities of daily living even while undergoing TNP therapy for a large diabetic foot wound.

Figure 16: Large plantar diabetic foot wound prior to Topical Negative Pressure dressing application

In a recent study,⁶⁷ we evaluated the potential ability of a removable cast walker (RCW) combined with a bridged TNP dressing technique to offload the plantar aspect of the foot. The combined device offloaded the foot nearly as well as the RCW alone. The difference between the TNP-RCW and the RCW amounted to only 2% of the initial barefoot pressure imparted to the wound (Figure 20). The results of this study suggested that a modified topical negative pressure dressing coupled with a robust RCW may not impart undue additional stress to the plantar foot. This may allow patients to retain some degree of protected freedom (and a potentially reduced length of hospital stay) while still allowing for the beneficial effects of TNP therapy and sufficient offloading.

Figure 17: Open-celled foam placed into large plantar diabetic foot wound

Figure 18: Clear topical negative pressure dressing applied. Note sponge bridge allowing for exit of dressing from dorsum of the foot. Tube is connected to V.A.C.® device

Figure 19: Topical negative pressure dressing with removable cast walker converted into "instant total contact cast", allowing for protected weight bearing. V.A.C.® device may be attached to waist-worn pack or backpack


Figure 20: Plantar pressure in bridged topical negative pressure device⁶⁷

RCW = removable cast walker

Bridged TNP = Topical Negative Pressure device applied plantarly with sponge dressing bridged to dorsum of foot

Another hybrid dressing has been developed to allow for more facility in treating heel wounds. This modality has the benefit of providing a simple, single heel cup which can completely surround the heel whilst allowing for bridging the tubing anteriorly (Figure 21). This provides the clinician the ability to offload the posterior pressure-related wound with the device of his or her choice.

Figure 21a: V.A.C.® GranuFoam™ (black, PU foam) is conformed to fit the heel

Figure 21b: TRAC™ pad bridge located on top of the foot for patient comfort

Conclusion

In conclusion, topical negative pressure therapy (TNP) appears to be a safe, efficacious, cost-effective, and practical method for promoting wound healing in persons with diabetes. Appropriate use of this modality requires careful patient selection and uncompromising attention to wound débridement, dressing application and concomitant pressure off-loading. We believe that this modality may facilitate our overall goal of improved rates of wound healing, and as a consequence, a reduction in lower extremity amputation.

Bibliography

- 1) Levin ME. Preventing amputation in the patient with diabetes. *Diabetes Care*. 1995; 18: 1383 – 1394.
- 2) Joslin EP. The menace of diabetic gangrene. *N. Engl. J. Med.* 1934; 211: 16 – 20.
- 3) Brand PW. The diabetic foot. In: Ellenberg M Rifkin H (ed.) *Diabetes Mellitus, Theory and Practice*. 3rd edn. New York: Medical Examination Publishing; 1983: 803 – 828.
- 4) Abbott CA, Vileikyte L, Williamson S, Carrington AL, Boulton AJM. Multicenter study of the incidence of and predictive risk factors for diabetic neuropathic foot ulceration. *Diabetes Care*. 1998; 21: 1071 – 5.
- 5) Abbott CA, Carrington AL, Ashe H, Bath S, et al. The North-West Diabetes Foot Care Study: incidence of, and risk factors for, new diabetic foot ulceration in a community-based patient cohort. *Diabetic Med.* 2002; 19: 377 – 84.
- 6) Boulton AJM, Vileikyte L. Diabetic foot problems and their management around the world. In: Bowker JH, Pfeifer MA (Eds) *Levin and O'Neal's The Diabetic Foot*. Mosby, St. Louis, 6th edn; 2001: 261 – 271.
- 7) Abbas ZG, Gill GV, Archibald LK. The epidemiology of diabetic limb sepsis: an African perspective. *Diabet Med.* 2002; 19: 895 – 9.
- 8) Tapp RJ, Shaw JE, de Courten MP, Dunstan DW, Welborn TA, Zimmet PA; AusDiab Study Group. *Diabetic Med.* 2003; 20: 105 – 13.
- 9) Young BA, Maynard C, Reiber G, Boyko EJ. Effects of ethnicity and nephropathy on lower-extremity amputation risk among diabetic veterans. *Diabetes Care*. 2003; 26: 495 – 501.
- 10) Boulton AJM. The pathway to ulceration. In: Boulton AJM, Connor H, Cavanagh PR (Eds) *The Foot in Diabetes*, Wiley, Chichester UK, 3rd Edn; 2000: 19 – 32.
- 11) Reiber GE, Vileikyte L, Boyko EJ, del Aguila M, Smith DG, Lavery LA, Boulton AJM. Causal pathways for incident lower-extremity ulcers in patients with diabetes from two settings. *Diabetes Care*. 1999; 22: 157 – 62.
- 12) Armstrong DG, Nguyen HC, Lavery LA, van Schie CH, Boulton AJ, Harkless LB. Off-loading the diabetic foot wound: a randomized clinical trial. *Diabetes Care*. 2001; 24: 1019 – 22.
- 13) Armstrong DG, Lavery LA, Kimbriel HR, Nixon BP, Boulton AJ. Activity patterns of patients with diabetic foot ulceration: Patients with active ulceration may not adhere to a standard pressure off-loading regimen. *Diabetes Care*. 2003; 26 (9): 2595 – 2597.
- 14) Boulton AJ, Armstrong DG. Trials in neuropathic diabetic foot ulceration: Time for a paradigm shift? *Diabetes Care*. 2003; 26 (9): 2689 – 2690.
- 15) American Diabetes Association. Consensus Development Conference on Diabetic Foot Wound Care. *Diabetes Care*. 1999; 22 (8): 1354.
- 16) Hogge J, Krasner D, Nguyen HC, Harkless LB, Armstrong DG. The potential benefits of advanced therapeutic modalities in the treatment of diabetic foot wounds. *J Am Podiatr Med Assoc*. 2000; 90 (2): 57 – 65.
- 17) Banwell PE. Topical negative pressure therapy in wound care. *J Wound Care*. 1999; 8: 79 – 84.
- 18) Banwell PE, Teot L. Topical negative pressure (TNP): the evolution of a novel wound therapy. *J Wound Care*. 2003; 12 (1): 22 – 28.
- 19) Argenta LC, Morykwas MJ. Vacuum-assisted closure: A new method for wound control and treatment: clinical experience. *Ann Plast Surg*. 1997; 38 (6): 563 – 576.
- 20) Morykwas MJ, Argenta LC, Shelton-Brown EI, McGuirt W. Vacuum-assisted closure: A new method for wound control and treatment: animal studies and basic foundation. *Ann Plast Surg*. 1997; 38 (6): 553 – 562.
- 21) Fleischmann W, Becker U, Bischoff M, Hoekstra H. Vacuum sealing: indication, technique and results. *Eur J Orthop Surg & Trauma*. 1995; 5: 37 – 40.
- 22) Banwell PE, Withey S, Holten IW. The use of negative pressure to promote healing. *Br J Plast Surg*. 1998; 51(1): 79.

- 23) Mullner T, Mrkonjic L, Kwasny O, Vecsei V. The use of negative pressure to promote the healing of tissue defects: a clinical trial using the vacuum sealing technique. *Br J Plast Surg*. 1997; 50: 194 - 199.
- 24) Banwell PE. Novel perspectives in wound care: topical negative pressure therapy. *European Tissue Repair Society Bulletin*. 2002; 9 (2): 49 - 50.
- 25) Gustafsson R, Johnsson P, Algotsson L, Blomquist S, Ingemansson R. Vacuum assisted closure therapy guided by C-reactive protein level in patients with deep sternal wound infection. *J Thorac Cardiovasc Surg*. 2002; 123 (5): 895 - 900.
- 26) Mooney JF3rd, Argenta LC, Marks MW, Morykwas MJ, DeFranzo AJ. Treatment of soft tissue defects in pediatric patients using the VAC system. *Clin Orthop*. 2000; 376: 26 - 31.
- 27) S M Jones, P E Banwell, P G Shakespeare, E Tiernan. Complications of topical negative pressure (TNP) therapy in small diameter wounds. *Plast Reconstr Surg*. 2003 (submitted).
- 28) Banwell PE, Evison D, Whitworth I. Vacuum therapy in degloving injuries of the foot: technical refinements. *Br J Plast Surg*. 2002; 553: 264 - 266.
- 29) Fleischman W, Russ M, Westhauser A. Vacuum sealing as a carrier system for controlled local drug administration in wound infection. *Unfallchirurg*. 1998; 101 (8): 649 - 54.
- 30) Bucalo B, Eaglstein WH, Falanga V. Inhibition of cell proliferation by chronic wound fluid. *Wound Rep Regen*. 1993; 1: 181 - 186.
- 31) Falanga V. Growth factors and chronic wounds: the need to understand the microenvironment. *J Dermatol*. 1992; 19: 667 - 72.
- 32) Wysocki AB, Grinnell F. Fibronectin profiles in normal and chronic wound fluid. *Lab Invest*. 1990; 63: 825 - 31.
- 33) Wysocki AB, Staiano-Coico L, Grinnell F. Wound fluid from chronic leg ulcers contains elevated levels of metalloproteinases MMP-2 and MMP-9. *J Invest Dermatol*. 1993; 101:64 - 8.
- 34) Reuler JB, Cooney TG. The pressure sore: pathophysiology and principles of management. *Ann Intern Med*. 1981; 94: 661 - 6.
- 35) Witkowski JA, Parish LC. Histopathology of the decubitus ulcer. *J Am Acad Dermatol*. 1982; 6:1014 - 21.
- 36) Ilizarov GA. The tension-stress effect on the genesis and growth of tissues: Part II. The influence of the rate and frequency of distraction. *Clin Orthop*. 1989: 263 - 85.
- 37) Thoma R. Ueber die histomechanik des gefasssystems und die pathogenese der angiosklerose. *Virchows Archiv F Path Anat*. 1911; 204: 1 - 74.
- 38) Olenius M, Dalsgaard CJ, Wickman M. Mitotic activity in expanded human skin. *Plast Reconstr Surg*. 1993; 91: 213 - 6.
- 39) Cherry GW, Austad E, Pasyk K, McClatchey K, Rohrich RJ. Increased survival and vascularity of random-pattern skin flaps elevated in controlled, expanded skin. *Plast Reconstr Surg*. 1983; 72: 680 - 7.
- 40) Ingber D. Integrins as mechanochemical transducers. *Curr Opin Cell Biol*. 1991; 3: 841 - 8.
- 41) Winter GD. Formation of the scab and the rate of epithelialization of superficial wounds in the skin of the young domestic pig. *Nature*. 1962; 193: 293.
- 42) Winter GD, Scales JT. Effect of air drying and dressings on the surface of a wound. *Nature*. 1963;197: 91.
- 43) Svensjo T, Pomahac B, Yao F, Slama J, Eriksson E. Accelerated healing of full-thickness skin wounds in a wet environment. *Plast Reconstr Surg*. 2000; 106: 602 - 12; discussion 613 - 4.
- 44) Field FK, Kerstein MD. Overview of wound healing in a moist environment. *Am J Surg*. 1994; 167:2S - 6S.
- 45) Breuing K, Eriksson E, Liu P, Miller DR. Healing of partial thickness porcine skin wounds in a liquid environment. *J Surg Res*. 1992; 52: 50 - 8.
- 46) Joseph E, Hamori CA, Bergman S et al. A prospective randomized trial of vacuum-assisted closure versus standard therapy of chronic nonhealing wounds. *Wounds* 2000;12 (3): 60 - 67.
- 47) Walgenbach KJ, Riabikhin AW, Andree C et al. Induction of angiogenesis following vacuum sealing. *ZfW V.A.C. Therapie* Nr.13/2 (2000).
- 48) Gouttefangeas C, Eberle M, Ruck P, et al. Functional T lymphocytes infiltrate implanted polyvinyl alcohol foams during surgical wound closure therapy. *Clin Exp Immunol*. 2001; 124: 398 - 405.
- 49) Simpson DM, Ross R. The neutrophilic leukocyte in wound repair a study with antineutrophil serum. *J Clin Invest*. 1972; 51: 2009 - 3.

- 50) Efron JE, Frankel HL, Lazarou SA, Wasserkrug HL, Barbul A. Wound healing and T-lymphocytes. *J Surg Res.* 1990; 48: 460 - 3.
- 51) Singer AJ, Clark RA. Cutaneous wound healing. *N Engl J Med.* 1999; 341: 738 - 46.
- 52) Thornton FJ, Schaffer MR, Barbul A. Wound healing in sepsis and trauma. *Shock.* 1997; 8: 391 - 401.
- 53) Fabian TS, Kaufman HJ, Lett ED, et al. The evaluation of subatmospheric pressure and hyperbaric oxygen in ischemic full-thickness wound healing. *Am Surg.* 2000; 66: 1136 - 43.
- 54) Philbeck TE, Jr., Whittington KT, Millsap MH, Briones RB, Wight DG, Schroeder WJ. The clinical and cost effectiveness of externally applied negative pressure wound therapy in the treatment of wounds in home healthcare Medicare patients. *Ostomy Wound Manage.* 1999; 45: 41 - 50.
- 55) Ferrell BA, Osterweil D, Christenson P. A randomized trial of low-air-loss beds for treatment of pressure ulcers. *Jama.* 1993; 269: 494 - 7.
- 56) Evans D, Land L. Topical negative pressure for treating chronic wounds: a systematic review. *Br J Plast Surg.* 2001; 54: 238 - 42.
- 57) Pham HT, Rosenblum GI, Lyons TD et al. Evaluation of a human skin equivalent for the treatment of diabetic foot ulcers in a prospective, randomized, clinical trial. *Wound.* 1999; 11 (4): 79 - 86.
- 58) McCallon SK, Knight CA, Valiulus JP, Cunningham MW, McCulloch JM, Farinas LP. Vacuum-assisted closure versus saline-moistened gauze in the healing of postoperative diabetic foot wounds. *Ostomy Wound Manage.* 2000; 46: 28 - 32, 34.
- 59) Scheufler O, Peek A, Kania NM et al. Problem-adapted application of vacuum occlusion dressings: Case report and clinical experience. *Eur J Plast Surg.* 2000; 23: 386 - 390.
- 60) Brand PW, Yancey P. *The Gift of Pain.* Grand Rapids: Zondervan, 1997.
- 61) Armstrong DG, Lavery LA, Harkless LB. Validation of a diabetic wound classification system. The contribution of depth, infection, and ischemia to risk of amputation see comments). *Diabetes Care.* 1998; 21 (5): 855 - 859.
- 62) Yuan-Innes MJ, Temple CL, Lacey MS. Vacuum-assisted wound closure: a new approach to spinal wounds with exposed hardware. *Spine.* 2001; 26 (3): E30 - 33.
- 63) Armstrong DG, Lavery LA, Abu-Rumman P, et al. Outcomes of subatmospheric pressure dressing therapy on wounds of the diabetic foot. *Ostomy Wound Manage.* 2002; 48 (4): 64 - 68.
- 64) Morykwas MJ, Argenta LC. Nonsurgical modalities to enhance healing and care of soft tissue wounds. *J South Orthop Assoc.* 1997; 6 (4): 279 - 288.
- 65) Greer SE, Duthie E, Cartolano B, Koehler KM, Maydick-Youngberg D, Longaker MT. Techniques for applying subatmospheric pressure dressing to wounds in difficult regions of anatomy. *J Wound Ostomy Continence Nurs.* 1999; 26 (5): 250 - 253.
- 66) Armstrong DG, Short B, Nixon BP, Boulton AJM. Technique for fabrication of an "instant" total contact cast for treatment of neuropathic diabetic foot ulcers. *J Amer Podiatr Med Assn.* 2002; 92: 405 - 408.
- 67) Armstrong DG, Kunze K, Short B, et al. Use of the instant total contact cast and a modified topical negative pressure (VAC) dressing technique to allow protected walking in persons with complex plantar diabetic foot wounds. *J Amer Podiatr Med Assn.* 2004 (submitted).

CURRENT TITLES

Advanced Management of the Open Abdomen

Skin Grafting: The Role of Topical Negative Pressure (TNP) Therapy

Topical Negative Pressure: Management of Complex Diabetic Foot Wounds

Surgical Wound Infection: The Role of Topical Negative Pressure Therapy

FORTHCOMING TITLES IN THIS SERIES

The Science of Topical Negative Pressure Therapy

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W. Lang¹
R. E. Horch²

Distal Extremity Reconstruction for Limb sSalvage in Diabetic Foot Ulcers with Pedal Bypass, Flap Plasty and Vacuum Therapy

Abstract

There are numerous factors resulting in tissue damage and foot ulcers in diabetic patients. Critical ischaemia should be routinely excluded and an examination of the limb perfusion must be done in all cases with clinical suspicion of reduced perfusion. Pedal bypass grafting procedures in combination with flap plasty techniques achieve good results in patients with non-healing diabetic ulcers. Timing of the procedure is essential to avoid progression of gangrene. Good results of pedal bypass grafts support an indication at an early stage. Vacuum assisted wound therapy and local wound surgery can be done first to prepare patients with septic foot lesions for two-staged procedures

Key words

Diabetic foot · pedal bypass · vacuum therapy · flap plasty · microvascular surgery · distal extremity reconstruction

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Consolidation of difficult wound areas at diabetic foot syndrome – a case report

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Summary

In the following I comment on a patient in advanced stage IV of paod with phlegmon of the right lower leg. Duration of hospital treatment was 98 days. After conservative therapy, new settings in the treatment of diabetes mellitus, hypertension and associated renal insufficiency a successful subintimal angioplasty (PTA) of the posterior tibial artery was carried out. Subsequent to that the vast local findings have been debrided as well as the 3rd and 4th toe have been resected. By stage adapted wound management as well as accompanying Prostavasin®-therapy and hyperbaric oxygen treatment consolidation of the extended wound could be achieved. Ongoing treatment under V.A.C.-therapy with the GranuFoam® results in good granulation process. The last two weeks of the vacuum therapy were continued outpatiently in arrangement with the health insurance and the KCI company in the local wound-centre Hannover. For covering the remained wound defect of the right foot autologous mesh graft plastic was carried out later on. Thereby the 2nd toe has been further resected. Postoperatively the mesh graft showed good ingrowing tendencies after two cycles of V.A.C.-therapy with the WhiteFoam®. After discharge consistently outpatient examinations confirmed these findings.

Key words

Vacuum-therapy, wound management, diabetic foot syndrome

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Schlüsselwörter: Diabetisches Fuss-Syndrom, V.A.C.® Therapy

12 years experience with V.A.C.® Therapy
in diabetic foot syndrome

Summary. Background: Negative pressure wound therapy (V.A.C.® Therapy) is proven as the most effective treatment after partial diabetic foot amputation.

Methods: During the last 12 years we treated 209 patients with diabetic foot lesions in stage II B, D and III B, D, respectively, according the *Armstrong-Score* in the SMZ Ost Donauspital, Vienna. Mean age was 64 (35–86) years and duration of V.A.C.® Therapy was 21 (15–34) days.

Results: In 95 patients (87.9%) major amputation was avoided using the negative pressure wound therapy.

Conclusions: In our hands V.A.C.® Therapy is proven to lower the number of major amputation in patients with diabetic foot syndrome.

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Keywords: Diabetic foot syndrome, V.A.C.® Therapy

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Comparison of Negative Pressure Wound Therapy Utilizing Vacuum-Assisted Closure to Advanced Moist Wound Therapy in the Treatment of Diabetic Foot Ulcers– A Multicenter Randomized Controlled Trial

Peter A. Blume, DPM¹

Jodi Walters, DPM²

Wyatt Payne, MD³

Jose Ayala, MD⁴

John Lantis, MD⁵

²Southern Arizona Veterans Affairs Medical Center, Surgery, Tucson, AZ

³Institute for Tissue Regeneration, Repair, and Rehabilitation, Bay Pines VA Healthcare System, Bay Pines, FL and The University of South Florida, Division of Plastic Surgery, Tampa, FL

⁴Valley Baptist Hospital, Brownsville, TX

⁵Saint Luke's Roosevelt, New York, NY

Running Title: VAC NPWT: Diabetic Foot Ulcer RCT

Corresponding Author:

Dr. Peter A. Blume

¹North American Center for Limb Preservation

506 Blake Street, New Haven, CT 06515

peter.b@snet.net

Received for publication 19 November 2007 and accepted in revised form 18 December 2007

ABSTRACT

Objective: To evaluate safety and clinical efficacy of Negative Pressure Wound Therapy (NPWT) compared to Advanced Moist Wound Therapy (AMWT) to treat diabetic patients with foot ulcers.

Research Design And Methods: This multicenter randomized controlled trial enrolled 342 patients mean age 58 years; 79% male. Complete ulcer closure was defined as skin closure (100% re-epithelization) without drainage or dressing requirements. Patients were randomized to either NPWT (Vacuum-Assisted Closure) or AMWT (predominately hydrogels and alginates) and received standard off-loading therapy as needed. The trial evaluated treatment until Day 112 or ulcer closure by any means. Patients whose wound achieved ulcer closure were followed at 3 and 9 months. Each study visit included closure assessment by wound exam and tracings.

Results: A greater proportion of foot ulcers achieved complete ulcer closure with NPWT (73/169, 43.2%) than AMWT (48/166, 28.9%) within the 112-day Active Treatment Phase ($p=0.007$). Kaplan-Meier median estimate for 100% ulcer closure was 96 days (95% CI: 75.0, 114.0) for NPWT and not determinable for AMWT ($p=0.001$). NPWT patients experienced significantly ($p=0.035$) fewer secondary amputations. The proportion of home care therapy days to total therapy days for NPWT was 9471/10579 (89.5%) and 12210/12810 (95.3%) for AMWT. In assessing safety, no significant difference between the groups was observed in treatment-related complications such as infection, cellulitis, and osteomyelitis at 6 months.

Conclusions: NPWT appears to be as safe as, and more efficacious, than AMWT for the treatment of diabetic foot ulcers.

This study was registered with ClinicalTrials.gov as NCT00432965.

In 2005, the Centers for Disease Control estimated the prevalence of Diabetes in the United States to be 20.8 million people.(1) A disabling complication with this disease is foot ulcer development (2, 3), which leads to non-healing chronic wounds that are difficult to treat. Moreover, diabetic foot ulcers (DFUs) are a significant risk factor for nontraumatic foot amputations in persons with diabetes.(4)

Various DFU treatments have been reported in the literature, including advanced moist wound therapy (AMWT) (5, 6), bioengineered tissue or skin substitutes (7, 8), growth factors (9, 10), electric stimulation (11), and negative pressure wound therapy (NPWT) (12). The treatment success is dependent on ulcer chronicity, patient compliance, appropriate off-loading of the appendage, and the therapy's mechanisms of action.

NPWT is a non-invasive system that creates a localized controlled subatmospheric (negative) pressure environment. In this study NPWT was provided by the Vacuum-Assisted Closure Therapy System (KCI USA, San Antonio, Texas) which promotes wound healing by delayed primary or secondary intention through creating a moist wound environment, preparing the wound bed for closure, reducing edema, and promoting granulation tissue formation and perfusion. Vacuum-Assisted Closure Therapy is indicated for use in all care settings and in patients with a variety of wound types including diabetic foot ulcers.

This multicenter randomized controlled trial (RCT) evaluated the safety and efficacy of NPWT compared to AMWT (predominately hydrogels and alginates) for the treatment of DFUs. The cost benefit analysis will be reported in a future publication.

RESEARCH DESIGN AND METHODS

The patient population consisted of diabetic adults ≥ 18 years with a Stage 2 or 3 (as defined by Wagner's Scale) calcaneal, dorsal or plantar foot ulcer $\geq 2 \text{ cm}^2$ in area after debridement.(13) Adequate blood circulation (perfusion) was assessed by Dorsum Transcutaneous Oxygen Test (TcPO_2) ≥ 30 mmHg, Ankle Brachial Index (ABI) ≥ 0.7 and ≤ 1.2 with toe pressure ≥ 30 mmHg, or Doppler arterial waveforms that were triphasic or biphasic at the ankle of the affected leg. Patients with recognized active Charcot disease; ulcers resulting from electrical, chemical or radiation burns; presence of untreated osteomyelitis, cellulitis, collagen vascular disease, or ulcer malignancy were excluded from the study. Patients with uncontrolled hyperglycemia ($\text{HbA}_{1c} > 12\%$) or inadequate lower extremity perfusion were not enrolled. Exclusion criteria also included ulcer treatment with normothermic or hyperbaric oxygen therapy; concomitant medications such as corticosteroids, immunosuppressive medications, or chemotherapy; recombinant or autologous growth factor products; skin and dermal substitutes within 30 days of study start; or use of any enzymatic debridement treatments. Pregnant or nursing mothers were excluded from study participation.

The NPWT system utilized in this study was Vacuum Assisted Closure Therapy. The system consists of 3 components: a negative pressure generating unit with a disposable canister, a pad with evacuation tube, and a reticulated, open cell sterile polyurethane or a dense open-pore polyvinyl alcohol foam dressing cut to fit the wound. The system unit is programmed to deliver controlled negative pressure ranging from 50 to 200 mmHg. NPWT was applied to the ulcer as specified by manufacturer's guidelines(14) and treated until ulcer closure, sufficient granulation tissue formation for healing by primary, secondary or surgical intention, or Day 112.

AMWT dressings were used according to Wound, Ostomy and Continence Nurses (WOCN) Society guidelines (6), and institutional treatment protocols, consistent with standards of care for treating DFUs. Skin substitutes, cytokines, recombinant human platelet-derived growth factors, or similar therapies as outlined in the exclusion criteria were not used in either group during the Active Treatment Phase (ATP).

The primary efficacy endpoint was incidence of complete ulcer closure. Secondary endpoints included reduction of ulcer surface area over time, time to achieve ulcer closure by either surgery or secondary intention, and reduction in complications, including secondary amputations. Complete ulcer closure was defined as skin closure (100% re-epithelization) without drainage or dressing requirements.

The sample size was based on a type I error probability set at 0.05, with 80% power. Detection of a 20% difference between treatment groups required 206 evaluable patients with a treatment-to-study ratio of 1:1. Sample size was set at 338 to account for subject withdrawal or loss to follow-up.

Randomization was accomplished by generating blocks of numbers through www.randomizer.org. Numbers were assigned to a treatment group and sealed in opaque envelopes containing black paper labeled with treatment and patient ID. Envelopes were sequentially numbered prior to clinical trial site distribution. At patient randomization, treatment was assigned based on the next sequentially labeled envelope.

FDA guidelines (15) state that in some devices, it is impractical or unethical to implement a control treatment that mimics the test product and allows masking. In this study, the physical differences between treatment regimens (eg, hydrogels and NPWT) can be so distinctive that it is not

possible to blind either the patient or physician to the treatment post randomization.

Study design. This study was a prospective RCT initiated at 37 diabetic foot and wound clinics, and hospitals. A total of 342 patients were enrolled at 1 Canadian and 28 USA sites. Prior to randomization, patients were screened for neuropathy, adequate perfusion, and glycemic control. All foot ulcers were assessed and debrided as needed within two days of randomization.

Patients were examined weekly for the first 4 weeks (Day 28) then every other week until Day 112, or ulcer closure by any means. At each study visit, ulcers were assessed for area via wound tracing, ulcer closure, and/or adequate granulation tissue formation. NPWT dressing changes were performed every 48 to 72 hours, no less than 3 times per week. AMWT randomized patients were treated based on manufacturer's guidelines. All patients received off-loading therapy as deemed necessary. Patients achieving ulcer closure were followed at 3 and 9 months.

Statistical analysis. The KCI Global Biometrics Group conducted the safety and effectiveness data analyses for this trial. Primary analysis was based on an intention-to-treat (ITT) dataset composed of all randomized patients who signed an informed consent and received at least one post-baseline treatment.

Continuous demographic variables (age, gender, etc.) were summarized for the study population as descriptive statistics (number, mean, Standard Deviation [SD]), median, minimum, and maximum values, and 95% two-sided confidence limits and compared between groups with a two-sample *t*-test or Wilcoxon Rank-Sum test.

An interim ITT analysis of the primary endpoint was conducted at the enrollment of 275 patients to confirm the consistency of trial design assumptions with observed data. The Lan-Demets group sequential boundaries (two-sided symmetric boundaries) were calculated

using the alpha-spending function of O'Brien-Fleming.

Data were analyzed based on differences in incidence of closure between treatment groups and assessed with a test of proportions, using chi-square or Fisher's Exact test, as appropriate. Treatment differences in ulcer area (cm^2) were calculated (ulcer area [cm^2] = ulcer length [cm] x ulcer width [cm] x $\pi/4$) and analyzed using Analysis of Covariance (ANCOVA) with baseline area used as a covariate. Median time to closure was based on number of days from baseline until closure via Kaplan-Meier (KM) survival analysis. Patients not achieving closure were censored using last date of observation.

RESULTS

Figure 1 describes the patient flow through each stage of this RCT including discontinued patients. During the course of the study, 384 patients were consented and screened for inclusion. Of these, 42 patients were excluded due to inclusion/exclusion criteria, patient refusal to participate, or withdrawal of consent, and 342 were enrolled. Seven patients did not receive treatment and 335 patients were analyzed. All patient data collected during ATP were included in the ITT analyses.

The data suggest that no statistically significant demographic differences existed between treatment arms (Table 1). The mean patient population age was 58 years and predominantly male (78.5%). Percentages of patients treated for study ulcer infections at baseline were 29.6% (50/169) for NPWT and 27.1% (45/166) for AMWT. The randomization method resulted in an even distribution of characteristics between treatment groups.

Patients receiving off-loading were 164/169 (97.0%) for NPWT and 162/166 (97.6%) for AMWT. Patients were treated in both acute and home care settings. The

proportion of home care therapy days to total therapy days was 9471/10579 (89.5%) for NPWT and 12210/12810 (95.3%) for AMWT.

Efficacy. Complete ulcer closure was defined as skin closure (100% re-epithelization) without drainage or dressing requirements. Within ATP, the NPWT group proportion was significantly ($p=0.007$) greater for complete ulcer closure than the AMWT group (73/169, 43.2% vs. 48/166, 28.9%). For patients completing ATP, completers analysis significantly ($p=0.001$) confirmed that a greater percentage of NPWT treated ulcers (60.8%, 73/120) achieved ulcer closure compared to AMWT (40.0%, 48/120). After sufficient wound bed preparation, 9.5% (16/169) NPWT-treated ulcers and 8.4% (14/166) AMWT-treated ulcers were surgically closed by split thickness skin grafts, flaps, sutures, or amputations. KM median time to complete ulcer closure was 96 days (95% CI: 75.0, 114.0) for NPWT ($p=0.001$). AMWT median time to complete ulcer closure could not be estimated (Figure 2). The NPWT mean duration of therapy was 63.6 days (SD 36.57) versus 78.1 days (SD 39.29) for AMWT.

To further evaluate NPWT effects, 75% ulcer closure, degree of granulation tissue formation, and ulcer area reduction were assessed. Significantly more NPWT patients (105/169, 62.1%) achieved 75% ulcer closure than AMWT patients (85/166, 51.2%; $p=0.044$). NPWT KM median estimate for 75% ulcer closure was 58 days (95% CI: 53.0, 78.0) and 84 days (95% CI: 58.0, 89.0) for AMWT ($p=0.014$). In assessing ulcer area, a significant difference between NPWT (-4.32 cm^2) and AMWT (-2.53 cm^2) from baseline was achieved on Day 28 ($p=0.021$).

The clinical treatment effect on wound bed preparation was assessed in 46 patients (24 NPWT, 22 AMWT), who presented with 0-10% granulation at baseline and achieved 76-100% granulation. Of these, 17/24 (70.8%) NPWT and 8/22 (36.4%) AMWT patients achieved 76-100% granulation tissue formation

($p=0.019$). NPWT KM median estimates for 76-100% granulation tissue formation were 56 days (95% CI: 42.0, 84.0) and 114 (95% CI: 44.0, -) for AMWT ($p=0.022$).

Safety. Table 2 reports treatment-related rates for secondary amputations, edema, wound infection, cellulitis, osteomyelitis, staphylococcal infection, and infected skin ulcers at 6 months. Significantly ($p=0.035$) fewer amputations were observed in NPWT patients (7/169, 4.1%) compared to AMWT patients (17/166, 10.2%). The majority of these amputations (2 and 13, respectively) were minor amputations. In all other categories, no significant differences were observed.

CONCLUSIONS

Results of the largest NPWT RCT to date demonstrate that NPWT is as safe as, and more efficacious, than AMWT in the treatment of DFUs. A significantly greater number of NPWT patients achieved complete ulcer closure and granulation tissue formation compared to AMWT patients. This was supported by a significant reduction in median time needed to heal DFUs. Approximately 90% of therapy days occurred in the home care setting for both treatments.

No significant difference was observed in ulcer related complications such as infection, cellulitis, and osteomyelitis. However, the study showed that AMWT patients had more than twice as many secondary amputations than those receiving NPWT.

Chronic DFUs present a significant challenge to treating physicians.(16) Treatment involves multiple modalities including debridement, assessment and treatment of infection, revascularization if indicated, and sufficient off-loading of the foot.(17) A key component of the healing process is debridement because it enables removal of devitalized and necrotic tissue.

Debridement is critically important to the initiation of healing. NPWT and other wound healing technologies work in conjunction with debridement as the foundation upon which the wound healing process can begin.(18) As observed in this clinical trial, the use of NPWT in concurrence with debridement of the affected foot increases the number of DFUs healed and decreases the length of time required for ulcer healing as compared to AMWT. In addition, the prescription of off-loading may have also contributed to positive results in both groups. Therefore, it appears that NPWT in addition to established standards of care enhances successful healing and closure of DFUs.

In this study, 14.3% more NPWT patients achieved complete ulcer closure in less median time to closure than those treated with AMWT. This parallels the findings by Argenta and Morykwas in 1997, in which they reported that the success of NPWT in chronic wounds is associated with removal of excess interstitial fluid, an increase in vascularity and associated decrease of bacterial colonization, and stimulation of granulation tissue formation through the response of wound tissue to the mechanical forces exerted by the application of negative pressure through the foam dressing.(10) In separate studies, Saxena (18) and Greene (19) have further elucidated the role of open pore foam dressing in the creation of micromechanical deformations of the wound surface. These micromechanical deformations are caused when negative pressure draws tissue into the foam pores. This stretches cells and promotes cell division that stimulates granulation tissue formation.(18)

As NPWT use has grown, the use of NPWT for wound bed preparation in conjunction with either delayed primary or secondary wound closure has also increased. In 2005, Armstrong and Lavery reported that the use of NPWT may be an alternative therapy to achieve an improved granulating wound bed in diabetic foot wounds in order to prepare the wound bed

for other closure techniques.(20) As shown in these data and previous studies, NPWT promotes granulation tissue formation, thereby allowing the clinician to determine course of closure.

DFUs are a significant risk for amputation.(4) In this study, the incidence of secondary amputations was significantly less ($p=0.035$) for NPWT (4.1%) compared to AMWT (10.2%). This finding substantiates other reports, in which diabetic foot wounds treated with NPWT trended towards fewer secondary amputations ($p=0.060$) than those treated with AMWT.(20) While the exact mechanism of the decrease in secondary amputations remains unclear, treatment of DFUs with NPWT appears to promote significant healing.

This study of 342 patients with DFUs showed that NPWT is as safe as and more efficacious than AMWT for the treatment of diabetic foot ulcers.

ACKNOWLEDGMENTS

KCI USA Incorporated (San Antonio, TX) supported this study. Conference abstracts reporting preliminary and interim study data were presented at the 2003 Symposium on Advanced Wound Care, 2004 World Union of Wound Healing Societies meeting, 2007 Symposium on Advanced Wound Care, and 2007 American Podiatric Medicine Association meeting. Investigators who enrolled patients for the V.A.C.[®] Therapy Diabetic Foot Ulcer Study were: Charles Anderson, MD, Madigan Army Hospital, Tacoma, WA; David Armstrong, DPM, North Chicago VA Medical Center and Rosalind Franklin University, North Chicago, IL; Jose Ayala, DPM, Valley Baptist Hospital, Brownsville,

TX; John Bennett, DPM, Des Moines University, Des Moines, IA; Scott Berman, MD, Tucson Vascular Surgery, Tucson, AZ; Peter Blume, DPM, North American Center for Limb Preservation, New Haven, CT; Milton Boden, MD, Aim Research, Atlanta, GA; William Bogey, MD, East Carolina University, Greenville, NC; Marc Dolce, DPM, Cleveland Foot and Ankle Clinic, Cleveland, OH; Kenneth Dolynchuk, MD, Saint Boniface General Hospital, Manitoba, Winnipeg, Canada; Tamara Fishman, DPM, Primary Foot Care Center, Incorporated, North Miami Beach, FL; Gabriel J. Halperin, DPM, Innovative Medical Technologies, Los Angeles, CA; Allen Holloway, MD, Maricopa Medical Center, Phoenix, AZ; John Lantis, MD, Saint Luke's Roosevelt, New York, NY; Brock Liden, DPM, Circleville Foot and Ankle, LLC, Circleville, OH; Jeffrey Page, DPM, Carl T. Hayden VA Medical Center, Podiatry Section, Phoenix, AZ; Marc Passman, MD, Vanderbilt University Clinical Trials Center, Nashville, TN; Wyatt Payne, MD, Bay Pines Veterans Affairs Medical Center, Bay Pines, FL; Alexander Reyzelman, DPM, Bay Area Foot Care, Castro Valley, CA; Bret Ribotsky, DPM, Podiatric Success, Inc., Boca Raton, FL; Kathleen Satterfield, DPM, Texas Diabetes Institute - University of Texas Health Science Center San Antonio, San Antonio, TX; Bhavesh Shah, DPM, South Texas Foot Institute, Poth, TX; David Skrobot, DPM, Genesis Health Care System, Zanesville, OH; Rodney Stuck, DPM, Hines Veterans Affairs Hospital, Hines, IL; Arthur Tallis, DPM, Hope Research Institute, Phoenix, AZ; Mike Vaardahl, DPM, Banner Health at North Colorado Medical Center, Greeley, CO; Jodi Walters, DPM, Southern Arizona Veterans Affairs Medical Center, Department of Surgery, Tucson, AZ; Robert Wunderlich, DPM, San Antonio, TX

REFERENCES

1. Centers for Disease Control and Prevention. National diabetes fact sheet: general information and national estimates on diabetes in the United States, 2005. Atlanta, GA: U.S. Department of Health and Human Services, Centers for Disease Control and Prevention, 2005
2. Ramsey SD, Newton K, Blough D, McCulloch DK, Sandhu N, Reiber GE, Wagner EH: Incidence, outcomes, and cost of foot ulcers in patients with diabetes. *Diabetes Care* 22: 382-387, 1999
3. Singh, N, Armstrong DG, Lipsky BA: Preventing foot ulcers in patients with diabetes. *JAMA* 293: 217-228, 2005
4. Brem H, Shehan P, Rosenberg HJ, Schnieder JS, Boulton AJM: Evidence-based protocol for diabetic foot ulcers. *Plast Reconstr Surg* 117: 193S-207S, 2006
5. Mulder G, Armstrong D, Seaman S: Standard, appropriate, and advanced care and medical-legal considerations: part one – diabetic foot ulcerations. *WOUNDS* 15: 92-106, 2003
6. Guideline for Management of Wounds in Patients with Lower-Extremity Arterial Disease. Glenview, IL: Wound Ostomy and Continence Nurses Society; 2002
7. Martson WA, Hanft J, Norwood P, Pollak R. The efficacy and safety of Dermagraft in improving the healing of chronic diabetic foot ulcers: results of prospective randomized trial. *Diabetes Care* 26: 1701-1705, 2003
8. Veves A, Falanga V, Armstrong DA, Sabolinski ML: Graftskin, a human skin equivalent, is effective in the management of noninfected neuropathic diabetic foot ulcers: a prospective randomized multicenter clinical trial. *Diabetes Care* 24: 2001-295, 2001
9. Robson MC, Payne WG, Garner WL, Biundo J, Giacalone VF, Cooper DM, Ouyang P: Integrating the results of phase IV (postmarketing) clinical trial with four previous trials reinforces the position that Regranex (becaplermin) gel 0.01% is an effective adjunct to the treatment of diabetic foot ulcers. *J Appl Res* 5: 35-45, 2005
10. Mannari RJ, Payne WG, Ochs DE, Walusimbi M, Blue M, Robson MC: Successful treatment of recalcitrant, diabetic heel ulcers with topical becaplermin (rhPDGF-BB) gel. *Wounds* 14: 116-121, 2002
11. Peters EJ, Lavery LA, Armstrong DG, Fleischli JG: Electric stimulation as an adjunct to heal diabetic foot ulcers: a randomized clinical trial. *Arch Phys Med Rehabil* 82:721-725, 2001
12. Argenta LC, Morykwas MJ: Vacuum-assisted closure: a new method for wound control and treatment: clinical experience. *Ann Plast Surg* 38:563-576, 1997
13. Wagner, FW: The dyvascular foot: a system for diagnosis and treatment. *Foot Ankle*, 2(2): 64-122, 1981
14. V.A.C.® Therapy Clinical Guidelines. A reference source for clinicians; KCI, San Antonio, TX; 2004
15. Guidance for Industry. Chronic Cutaneous Ulcer and Burn Wounds — Developing Products for Treatment (Draft Guidance). Rockville, MD: U.S. Department of Health and Human Services; 2000. 15
16. Sibbald RG, Torrance G, Hux M, Attard C, Milcovich N: Cost-effectiveness of becaplermin for non-healing neuropathic diabetic foot ulcers. *Ostomy Wound Manage* 43: 76-84, 2003
17. Hopf HW, Humphrey LM, Puzifferri N, West JM, Attinger CE, Hunt TK: Adjuncts to preparing wounds for closure: hyperbaric oxygen, growth factors, skin substitutes, negative pressure wound therapy (vacuum assisted closure). *Foot Ankle Clin* 6: 661-682, 2001

18. Saxena V, Hwang CW, Huang S, Eichbaum Q, Ingber D, Orgill DP: Vacuum-assisted closure: Microdeformations of wounds and cell proliferation. *Plast Reconstr Surg* 114: 1086-1096, 2004
19. Greene AK, Puder M, Roy R, Arsenault D, Kwei S, Moses MA, Orgill DP: Microdeformational wound therapy: Effects on angiogenesis and matrix metalloproteinases in chronic wounds of 3 debilitated patients. *Ann Plast Surg* 56: 418-422, 2006
20. Armstrong DF, Lavery LA. Negative pressure wound therapy after partial diabetic foot amputation: a multicentre, randomised controlled trial. *Lancet* 366: 1704-1710, 2005

TABLE 1. Patient demographics

Characteristics	NPWT* n=169	AMWT n=166
Age (years)	58 ± 12	59 ± 12
Sex (M/F)	141/28 (83/17)	122/44 (73/27)
Race		
African American	28 (16.6)	22 (13.3)
Caucasian	95 (56.2)	100 (60.2)
Hispanic	41 (24.3)	40 (24.1)
Native American	3 (1.8)	3 (1.8)
Other	2 (1.2)	1 (0.6)
Weight (kg)	99.2 ± 25.1	93.8 ± 25.6
Height (cm)	175.0 ± 9.6	175.0 ± 12.4
Current Smoker	34 (20.1)	32 (19.4)
Currently Use Alcohol	37 (21.9)	45 (27.1)
Type of Diabetes		
Type 1	15 (8.9)	14 (8.4)
Type 2	154 (91.1)	152 (91.6)
Pre-albumin (g/L)	21.1 ± 7.6	19.9 ± 7.9
Albumin (g/L)	3.4 ± 0.6	3.4 ± 0.8

Characteristics	NPWT* n=169	AMWT n=166
HbA _{1c}	8.3 ± 2.0	8.1 ± 1.9
ABI (mmHg)	1.0 ± 0.2	1.0 ± 0.2
Transcutaneous oxygen tension (TcPO ₂ ; mmHg)	43.2 ± 10.4	43.3 ± 12.5
Loss of Protective Sensation [†]	150 (90.4)	143 (88.8)
Ulcer Duration Prior to Treatment (days)	198.3 ± 323.5	206.0 ± 365.9
Baseline Wound Area (cm ²)	13.5 ± 18.2	11.0 ± 12.7
Received Off-Loading Therapy	164 (97.0)	162 (97.6)
Treated for Ulcer Infection Prior to Randomization	50 (29.6)	45 (27.1)
Therapy Received at Treatment Initiation		
NPWT	169 (100)	
Hydrogel		78 (47.0)
Alginate		31 (18.7)
Other		28 (16.9)
Saline		17 (10.2)
Collagen		11 (6.6)
Hydrocolloid		1 (0.6)

*Data are means ± SD or n (%) [†]Percentage based on available data

TABLE 2. Results of safety analysis

MeDRA System Organ Class	NPWT n=169 (%)	AMWT n=166 (%)	P-Value*
Secondary Amputations	7 (4.1)	17 (10.2)	0.035
Edema	5 (3.0)	7 (4.2)	0.571
Wound Infection	4 [†] (2.4)	1 [‡] (0.6)	0.371
Cellulitis	4 (2.4)	1 (0.6)	0.371
Osteomyelitis	1 (0.6)	0 (0.0)	--
Staphylococcus Infection	1 (0.6)	0 (0.0)	--
Infected Skin Ulcer	1 (0.6)	2 (1.2)	0.620

*Fisher's Exact Test

[†]1 moderate, 3 mild

[‡]1 moderate

FIGURE 1



FIGURE 2

Outcomes of Subatmospheric Pressure Dressing Therapy on Wounds of the Diabetic Foot

David G. Armstrong, DPM
Lawrence A. Lavery, DPM, MPH
Patricia Abu-Rumman, DPM
Eric H. Espensen, DPM
Jefferey R. Vazquez, DPM
Brent P. Nixon, DPM
Andrew J.M. Boulton, MD

Ostomy/Wound Management
April 2002;48(4):64-68

Reprinted with permission from *Ostomy/Wound Management*
Copyright © 2001, HMP Communications

This article was written by an independent clinician not employed by KCI. While the article supports general wound healing indications, it may also support indications that have been cleared by the FDA. Consult a physician or product labeling for proper indications, contraindications, and precautions for V.A.C.® therapy.

OUTCOMES OF SUBATMOSPHERIC PRESSURE DRESSING THERAPY ON WOUNDS OF THE DIABETIC FOOT

— David G. Armstrong, DPM; Lawrence A. Lavery, DPM, MPH; Patricia Abu-Rumman, DPM; Eric H. Espensen, DPM; Jefferey R. Vazquez, DPM; Brent P. Nixon, DPM; Andrew J.M. Boulton, MD

The purpose of this retrospective study was to evaluate outcomes of people with large diabetic foot wounds treated with subatmospheric pressure dressing therapy immediately following surgical wound debridement. Data were abstracted from the medical records of 31 consecutive patients with diabetes, 77.4% male (n = 24), aged 56.1 ± 11.7 years, presenting for care at two large multidisciplinary wound care centers. All patients received surgical debridement for indolent diabetic foot wounds and were subsequently started on a regimen of subatmospheric pressure dressing therapy delivered using a vacuum-assisted closure device for a mean of 4.7 ± 4.2 weeks (mode = 2 weeks) using a protocol that called for cessation of therapy when the wound bed approached 100% coverage with granulation tissue with no exposed tendon, joint capsule, or bone. Outcomes evaluated included time to complete wound closure, proportion of patients achieving wound healing at the level of initial debridement, and complications associated with use of the device. The mean duration of wounds before therapy was 25.4 ± 23.8 weeks. In patients treated with subatmospheric pressure dressing therapy, 90.3% (n = 28) of wounds healed

at the level of debridement without the need for further bony resection in a mean 8.1 ± 5.5 weeks. The remaining 9.7% (n = 3) went on to higher level amputation (below knee amputation = 3.2%, [n = 1] and transmetatarsal amputation = 6.5% [n = 2]). Complications included periwound maceration (19.4% [n = 6]), periwound cellulitis (3.2% [n = 1]), and deep space infection (3.2% [n = 1]). The authors concluded that appropriate use of subatmospheric pressure dressing therapy to achieve a rapid granular bed in diabetic foot wounds may have promise in treatment of this population at high risk for amputation and that a large, randomized trial is now indicated.

Ostomy/Wound Management 2002;48(4):64-68

Diabetic foot wounds commonly proceed to lower extremity amputation. They have been identified as the key component in nearly nine in 10 of these ablative procedures.¹ Clearly, healing these wounds in as rapidly a manner as possible may result in a substantial reduction in morbidity.

Dr. Armstrong is Director of Research in the Department of Surgery Southern Arizona Veterans Affairs Medical Center, Tucson, Ariz., and Visiting Senior Lecturer in the Department of Medicine, The Manchester Royal Infirmary, United Kingdom. Dr. Lavery is Associate Professor, Department of Orthopaedics, Loyola University, Chicago, Ill. Dr. Abu-Rumman, Dr. Espensen, and Dr. Vazquez are staff in the Department of Surgery, Southern Arizona Veterans Affairs Health Center. Dr. Boulton is Professor, University Department of Medicine, The Manchester Royal Infirmary, United Kingdom. Please address correspondence to: DG Armstrong, Department of Surgery, 3601 South Sixth Avenue, Tucson, AZ 85723; email: armstrong@usa.net.

One recently introduced method used to heal defects of the skin involves topical subatmospheric pressure dressing (SPD) therapy. This modality is most commonly delivered to the wound by way of an open-cell foam dressing covered with a clear membrane over the wound. The dressing is attached to a pump that provides either intermittent or continuous suction with a standard setting of 125 mm Hg.²⁴ Several studies have reported varying levels of success in different categories of wounds in both human and animal models.⁹⁻²² Although one randomized trial of 10 patients with diabetic foot wound exists,⁵ the authors are unaware of any other studies in the peer-reviewed medical literature that have specifically assessed outcomes (by way of primary data) associated with this type of therapy in this very high risk population. Therefore, this preliminary study was conducted to evaluate outcomes associated with application of SPD therapy on diabetic foot wounds.

Materials and Methods

Data were collected from 31 consecutive patients with diabetes (see Table 1), 77.4% male, aged 56.1 ± 11.7 years presenting for care at two large, multidisciplinary, referral-based wound care centers. All patients had a previous diagnosis of diabetes mellitus as assessed by their attending physician. All patients received surgical debridement for indolent diabetic foot wounds and were subsequently placed on SPD therapy delivered using VAC devices (Vacuum Assisted Closure, KCI, Inc., San Antonio, Texas, USA). Wounds tended to be large (27.9 ± 19.5 cm²), as they were evaluated after the aforementioned surgical debridement. Data were collected from patients over a 1-year period. Retrospective collection of these data was determined to be exempt from consent policy according to the governing institutional review board.

TABLE 1
DESCRIPTIVE CHARACTERISTICS OF PATIENTS

Number of patients	31
Age (years)	56.1 ± 11.7
% Male	77.4 (n = 24)
Wound duration before therapy (weeks)	25.4 ± 23.8
Surface area (cm ²)	27.9 ± 19.5
Depth (% University of Texas grade 1, 2, 3 respectively)	3.2 (n = 1), 45.2 (n = 14), 51.6 (n = 16)

University of Texas Grades

0. Pre or post-ulcerative wound, completely epithelialized
1. Superficial, not including tendon, capsule or bone
2. Exposed tendon, and/or closed joint capsule
3. Bone, or open joint exposed

All wounds were described and classified using the University of Texas diabetic foot wound classification system.²⁵ The prevalence of University of Texas grade 1, 2, and 3 lesions in this population was 3.2%, 45.2%, and 51.6%, respectively. The foot wounds were located in the digit/ray, transmetatarsal, midfoot, and heel (see Table 2). Although all wounds were initially infected, all were converted to stage "A" (noninfected, nonischemic) before device application, using a combination of surgical debridement and culture-directed systemic antibiotic therapy. Outcomes evaluated included time to complete wound closure, proportion of patients achieving wound healing at the level of initial debridement, and complications associated with use of the device. The mean duration of wounds before therapy was 25.4 ± 23.8 weeks. Following sharp debridement, patients were treated with SPD therapy at 125 mm Hg using a

Ostomy, Wound Management 2002;48:41-64 58

KEY POINTS

- To evaluate outcomes of subatmospheric pressure therapy in the clinical management of foot ulcers in people with diabetes, a chart review of 31 patients from two diabetic foot clinics was conducted.
- After an average of 4.7 weeks, most wounds were found to contain enough healthy granulation tissue to discontinue treatment. Periwound maceration was the most commonly observed side effect.
- Prospective, controlled clinical studies are needed to ascertain optimal treatment times and to obtain effectiveness and cost-effectiveness data.

TABLE 2
LOCATION OF WOUNDS

Location	%	N
Digit/Ray	51.6	16
Transmetatarsal	12.9	4
Midfoot	12.9	4
Heel	22.6	7
Total	100.0	31

protocol that called for consideration for cessation of therapy when the consistency of the wound bed approached 100% granulation tissue

with no overtly exposed tendon, joint capsule, or bone. Subatmospheric pressure dressing therapy dressings were changed every 48 hours.

All descriptive data are described as mean \pm standard deviation. The authors used Pearson's test to evaluate correlation between continuous variables such as time to healing based on wound duration. To compare differences in time to healing by gender, a Student's *t*-test for independent samples was used. For all analyses, a significance (α) level of 0.05 was used.²⁴

Results

Subatmospheric pressure dressing therapy was used for a mean 4.7 ± 4.2 weeks (mode = 2 weeks) until the wound bed approached 100% granular tissue (see Figure 1). No significant difference in age, gender, duration of wounds, or time to healing between centers was noted. In total, 90.3% ($n = 28$) of wounds healed at the level of debridement without the need for further bony resection in a mean 8.1 ± 5.5 weeks. The remaining 9.7% ($n = 3$) went on to higher level amputation (below knee amputation = 3.2% [$n = 1$] and transmetatarsal amputation = 6.5% [$n = 2$]). Complications included periwound maceration (19.4% [$n = 6$]), periwound cellulitis (3.2% [$n = 1$]), and deep space infection (3.2% [$n = 1$]). In this small sample, no significant difference in healing time ($P = 0.2$) or proportion of patients healing ($P = 0.6$) based on gender was noted. Likewise, no significant association was reported between age ($P = 0.4$) or duration of the wound ($P = 0.9$) and time to wound healing ($P = 0.4$). Similarly, the authors could detect no association between wound size and time to healing ($P = 0.5$).

Discussion

The use of SPD therapy has been proposed as a

novel method of manipulating the chronic wound environment in a way that potentially reduces bacterial burden and chronic interstitial wound fluid, increases vascularity and cytokine expression, and, to an extent, mechanically exploits the viscoelasticity of periwound tissues.^{24,25-28} Although SPD therapy has appeared anecdotally promising, only one small randomized trial evaluating it in treatment of diabetic foot wounds can be found in the literature. McCallon et al³ enrolled 10 patients with diabetes and concomitant foot ulceration. Patients in this study were evenly randomized (by way of an alternating randomization technique) to either SPD therapy or to a control group (twice daily dressing changes with saline-moistened gauze). At the conclusion of the trial, the authors reported a difference in time to healing in patients receiving SPD therapy compared with standard therapy (23 ± 17 vs. 43 ± 32 days). No statistical analysis was performed on this admittedly underpowered sample. Most of the patients in the SPD group received definitive closure by way of delayed primary closure (4 out of 5 SPD versus 2 out of 5 control). Although this study may have had identifiable issues related to selection and performance bias between groups (as no readily identifiable method of randomization occurred), it did suggest a guardedly optimistic assessment of this treatment regime. Further works in this area may help to substantiate the potential of SPD as noted in the McCallon study.

In the authors' study, which intended to provide pilot data for future randomized controlled trials, valuable experience was gained in determining how


Figure 1
Duration of subatmospheric pressure therapy (weeks).
The mean time of therapy was 4.7 ± 4.2 weeks
(mode = 2 weeks).

this technology might be best used in the diabetic foot. It has become the authors' general protocol to treat patients with SPD therapy until they reach coverage of vital periarticular, tendinous, or osseous structures and/or a 100% granular base. The clinician may then choose to employ other wound healing modalities to move the wound toward complete epithelialization. This is reflected in the relatively short amount of time the device was used in this study. The protocol used most frequently called for SPD therapy for 1 to 3 weeks (see Figure 1) to achieve the aforementioned goals. It also should be noted that the patients in this study had proportionately much larger wounds than those in the McCallon trial, largely because the centers participating in the study utilized SPD therapy following a protocol of aggressive surgical debridement of chronic or infected diabetic foot wounds.

Limitations

This was a retrospective study and, as such, could not control for the myriad variables that could confound many conclusions. For example, the authors were unable to control for vascular status, degree of glucose control, nutritional status, activity level, or numerous other potential variables that could be important in assessing overall results. However, the objective of this study was to report outcomes and experiences using this treatment at busy diabetic foot

clinics, not to provide evidence to show SPD is more (or less) effective than other modalities.

Conclusion

This study provides additional pilot data to support, in a guarded sense, the judicious use of subatmospheric pressure dressing therapy in the diabetic foot. The results of this study suggest that appropriate use of subatmospheric pressure dressing therapy to achieve a rapid granular base in diabetic foot wounds may show promise in treating this population of patients who are at high risk for lower extremity amputation. To the authors' knowledge, this is the largest study to date to evaluate vacuum therapy in this specific wound category. The study will provide pilot data for future large, multicenter, randomized controlled trials of this wound healing modality to further clarify its utility in treating indolent diabetic foot wounds. Ultimately, this may yield improved results in therapy and a reduction in the unnecessarily high prevalence of lower extremity amputations in the developed and developing world. - OWM

References

1. Pecoraro RE, Reiber GE, Burgess EM. Pathways to diabetic limb amputation: basis for prevention. *Diabetes Care*. 1990;13:513-521.
2. Argenta LC, Morykwas MJ. Vacuum-assisted closure: a new method for wound control and treatment: clinical experience. *Ann Plast Surg*. 1997;38(6):563-576; discussion 577.

3. Bowler PG, Duerden BI, Armstrong DG. Wound microbiology and associated approaches to wound management. *Clin Microbiol Rev.* 2001;14(2):244-269.
4. Morykwas MJ, Argenta LC. Nonsurgical modalities to enhance healing and care of soft tissue wounds. *J South Orthop Assoc.* 1997;6(4):279-288.
5. McCallon SK, Knight CA, Valiulus JP, Cunningham MW, McCulloch JM, Farinas LP. Vacuum-assisted closure versus saline-moistened gauze in the healing of postoperative diabetic foot wounds. *Ostomy/Wound Management.* 2000;46(8):28-34.
6. Greer SE, Longaker MT, Margiotta M, Mathews AJ, Kasabian A. The use of subatmospheric pressure dressing for the coverage of radial forearm free flap donor-site exposed tendon complications. *Ann Plast Surg.* 1999;43(5):551-554.
7. Morykwas MJ, Argenta LC, Shelton-Brown EI, McGuirt W. Vacuum-assisted closure: a new method for wound control and treatment: animal studies and basic foundation. *Ann Plast Surg.* 1997;38(6):553-562.
8. Greer SE, Duthie E, Cartolano B, Koehler KM, Maydick-Youngberg D, Longaker MT. Techniques for applying subatmospheric pressure dressing to wounds in difficult regions of anatomy. *Journal of Wound Ostomy Continence Nursing.* 1999;26(5):250-253.
9. Philbeck TE, Jr., Whittington KT, Millsap MH, Briones RB, Wight DG, Schroeder WJ. The clinical and cost effectiveness of externally applied vacuum wound therapy in the treatment of wounds in home healthcare Medicare patients. *Ostomy/Wound Management.* 1999;45(11):41-50.
10. Deva AK, Siu C, Nettle WJ. Vacuum-assisted closure of a sacral pressure sore. *Journal of Wound Care.* 1997;6(7):311-312.
11. Fabian TS, Kaufman HJ, Lett ED, et al. The evaluation of subatmospheric pressure and hyperbaric oxygen in ischemic full-thickness wound healing. *Am Surg.* 2000;66(12):1136-1143.
12. Banwell P, Withey S, Holten I. The use of vacuum to promote healing. *Br J Plast Surg.* 1998;51(1):79.
13. Morykwas MJ, David LR, Schneider AM, et al. Use of subatmospheric pressure to prevent progression of partial-thickness burns in a swine model. *J Burn Care Rehabil.* 1999;20(1 Pt 1):15-21.
14. Genecov DG, Schneider AM, Morykwas MJ, Parker D, White WL, Argenta LC. A controlled subatmospheric pressure dressing increases the rate of skin graft donor site reepithelialization. *Ann Plast Surg.* 1998;40(3):219-225.
15. Mullner T, Mrkonjic L, Kwasny O, Vecsei V. The use of vacuum to promote the healing of tissue defects: a clinical trial using the vacuum sealing technique. *Br J Plast Surg.* 1997;50(3):194-199.
16. Deva AK, Buckland GH, Fisher E, et al. Topical vacuum in wound management. *Med J Aust.* 2000;173(3):128-131.
17. Josty IC, Ramaswamy R, Laing JH. Vacuum assisted closure: an alternative strategy in the management of degloving injuries of the foot. *Br J Plast Surg.* 2001;54(4):363-365.
18. Tang AT, Ohri SK, Haw MP. Novel application of vacuum assisted closure technique to the treatment of sternotomy wound infection. *Eur J Cardiothorac Surg.* 2000;17(4):482-484.
19. Yuan-Innes MJ, Temple CL, Lacey MS. Vacuum-assisted wound closure: a new approach to spinal wounds with exposed hardware. *Spine.* 2001;26(3):E30-33.
20. Mooney JE, 3rd, Argenta LC, Marks MW, Morykwas MJ, DeFranzo AJ. Treatment of soft tissue defects in pediatric patients using the V.A.C. system. *Clin Orthop.* 2000(376):26-31.
21. Greer SE, Adelman M, Kasabian A, Galiano RD, Scott R, Longaker MT. The use of subatmospheric pressure dressing therapy to close lymphocutaneous fistulas of the groin. *Br J Plast Surg.* 2000;53(6):484-487.
22. Greer SE, Grossi EA, Chin D, Longaker MT. Subatmospheric pressure dressing for saphenous vein donor-site complications. *Ann Thorac Surg.* 2001;71(3):1038-1040.
23. Armstrong DG, Lavery LA, Harkless LB. Validation of a diabetic wound classification system: the contribution of depth, infection, and vascular disease to the risk of amputation. *Diabetes Care.* 1998;21:855-859.
24. Kirkwood BR. *Essentials of Medical Statistics.* Oxford, UK: Blackwell; 1988.
25. Fleischmann W, Lang E, Kinzl L. Vacuum assisted wound closure after dermatofasciotomy of the lower extremity. *Unfallchirurg.* 1996;99(4):283-287.
26. Sumpio BE, Banes AJ, Link GW, Iba T. Modulation of endothelial cell phenotype by cyclic stretch: inhibition of collagen production. *J Surg Res.* 1990;48(5):415-420.
27. Sumpio BE, Banes AJ, Link WG, Johnson G, Jr. Enhanced collagen production by smooth muscle cells during repetitive mechanical stretching. *Arch Surg.* 1988;123(10):1233-1236.
28. Ingber DE. Control of capillary growth and differentiation by extracellular matrix. Use of a tensegrity (tensional integrity) mechanism for signal processing. *Chest.* 1991;99(3 Suppl):34S-40S.