

V.A.C. ULTA™ NEGATIVE PRESSURE WOUND THERAPY WITH V.A.C. VERAFLU™ THERAPY



Topical Wound Solutions

The following topical solutions are compatible with V.A.C. VERAFLU™ Therapy based on KCI in-house testing of the V.A.C. VERAFLU™ Therapy disposables with respect to mechanical properties, biocompatibility, and solution interaction.*



Generic Solution Class	Trade Name	Considerations for Use with V.A.C. VERAFLU™ Therapy
Hypochlorite-based solutions (e.g. hypochlorous acid, sodium hypochlorite)	Dakin's Solution (quarter strength), Dermacyn®, Microcyn®	<ul style="list-style-type: none"> Based on results of compatibility testing, Dakin's Solution should not be used in concentrations greater than 0.125% (quarter strength). Consider using the fewest irrigation cycles and minimizing hold times to the lowest level that is clinically relevant.
Silver nitrate (0.5%)	Various	Silver nitrate is light sensitive. Protect V.A.C. VERALINK™ Instillation Tubing from light during use of silver nitrate.
Sulfur-based solutions (sulfonamides)	Mafenide acetate, Sulfamylon®	Reconstitute powder per manufacturer's labeling.

*KCI is solely responsible for the solution compatibility information provided herein. Always refer to the solution manufacturer's prescribing information when deciding whether to administer any solution with the V.A.C. ULTA™ Therapy System and for appropriate safety, efficacy and dosing information. Listing of specific solutions is not an endorsement of a specific solution or an indication of clinical efficacy or safety.



For more information on the benefits of V.A.C. VERAFLU™ Therapy, contact your **KCI Representative**.

NOTE: Specific indications, contraindications, warnings, precautions and safety information exist for KCI products and therapies. Please consult a physician and product instructions for use prior to application. Rx only.

V.A.C.ULTA™ Therapy System is indicated for use with solutions labeled for topical use.

©2011-2016 KCI Licensing, Inc. All rights reserved. All trademarks designated herein are proprietary to KCI Licensing, Inc. DSL#16-0362.US • LIT#29-A-212 • (Rev. 4/16)



Topical Wound Solutions (continued)

The following topical solutions are compatible with V.A.C. VERAFLU™ Therapy based on KCI in-house testing of the V.A.C. VERAFLU™ Therapy disposables with respect to mechanical properties, biocompatibility, and solution interaction.*



Generic Solution Class	Trade Name	Considerations for Use with V.A.C. VERAFLU™ Therapy
Biguanides (Polyhexanide)	PRONTOSAN® Wound Irrigation Solution	<p>1,000ml bottle with spike adaptor available through KCI (Part Number: 400230)</p> <p>Other solution container configurations can be:</p> <ul style="list-style-type: none"> spiked directly connected using a spike adapter if available transferred to a container having a spike port
Lidocaine HCl	Various, e.g. Xylocaine 2%	<ul style="list-style-type: none"> For compatibility assessment, lidocaine HCl was tested in final concentrations equal to or less than 0.1% in saline. Toxicity concerns may exist; consult with your pharmacist and/or prescribing physician for solution and patient-specific considerations when using lidocaine HCl as an additive to saline instillation solutions. No device-related considerations for use with V.A.C. VERAFLU™ Therapy
Isotonic solutions	Normal Saline for Irrigation, USP Lactated Ringers Irrigation	May need to be connected using a spike adapter.
Acetic acid	Acetic Acid Irrigation USP	For compatibility assessment, acetic acid was tested in final concentrations equal to or less than 0.25% in sterile water

*KCI is solely responsible for the solution compatibility information provided herein. Always refer to the solution manufacturer's prescribing information when deciding whether to administer any solution with the V.A.C.ULTA™ Therapy System and for appropriate safety, efficacy and dosing information. Listing of specific solutions is not an endorsement of a specific solution or an indication of clinical efficacy or safety.