

ARJO HUNTLEIGH
GETINGE GROUP

MINUET®
Nursing and Community Care Bed



...with people in mind

PRODUCT SPECIFICATIONS

The MINUET low-height nursing and community care bed provides a real solution for independent living.

In addition to the 230mm low-height, the MINUET bed comes with easy-to-use electric controls. The unique design features also benefit the carer by reducing the amount of physical positional adjustment required for the user.

Standard features

- Height range of 230mm - 630mm (100mm castors fitted)
- Electric backrest and height adjustment (2-section model)
- Electric backrest, height adjustment and knee-break (4-section model)
- Trendelenburg / reverse trendelenburg tilt
- 180kg safe working load
- Mechanical calf adjustment
- Removable sheet metal platform
- Easy to clean ABS moulded head and foot panels
- Ivory coloured frame

Accessories

- Clamp-on 3/4 length folding safety sides
- Self-help pole with handle (stainless steel)
- Range of pressure relieving and reducing mattresses



230mm low-height

Technical details

Safe working load	180 kg	
Overall length	2200 mm	
Overall width	920 mm	
Height variation	230 mm - 630 mm	
Backrest angle	0 - 70°	
Leg elevation	(2-section)	to 15°
	(4-section)	to 18°
Tilt (tilt model only)	head down	10°
	head up	10°
Product weight	79 kg	

Electrical data

Power input	1A max at 230V~,50Hz/60Hz	
Power output	24V DC	
Electrical shock protection	Class II Type B	
EMC	Complies with BC EN 60601-1-2 2001	

Optional features

- Electric knee-break adjustment
- Range of panel infill colours



trendelenburg/reverse trendelenburg tilt

ARJOHUNTLEIGH
GETINGE GROUP

ArjoHuntleigh Australia
PO Box 330
Hamilton Hill, Western Australia 6963
Phone: +61 (0) 8 9337 4111
Fax: +61 (0) 8 9337 9077
www.ArjoHuntleigh.com

...with people in mind

Only ArjoHuntleigh designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by ArjoHuntleigh. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice.

Minuet is a trademark belonging to the ArjoHuntleigh group.
© ArjoHuntleigh, 2008

MEMBER OF THE GETINGE GROUP