

PALAMED®

Medium viscosity bone cement. From the PALACOS® family.


TRUST IN GREEN

PALAMED® is a medium viscosity bone cement based on the proven PALACOS®. It contains the same high quality raw materials as the Gold Standard PALACOS® and has very good mechanical properties. Characteristically green in colour, PALAMED® contrasts optically with the surrounding tissue and thus facilitates easier handling during surgery. PALAMED® is ideal for use with vacuum mixing systems and suitable for large and medium-sized joints.

BENEFITS OF ANTIBIOTIC PROPHYLAXIS

PALAMED® bone cement is also available with added gentamicin for antibiotic prophylaxis under the name PALAMED® G. The local release of the antibiotic minimizes the risk of infection in the implant with reduced systemic load¹.

SPECTRUM OF ACTIVITY OF GENTAMICIN



ADVANTAGES AT A GLANCE


- Suitable for large and medium-sized joints
- Ideal for vacuum mixing
- Antibiotic prophylaxis from added Gentamicin
- Consistent quality based on proven PALACOS® raw materials
- Available with or without Gentamicin



¹ Engesaeter, L.B. et al.: Antibiotic prophylaxis in total hip arthroplasty, University of Bergen, Norway 2003.

PALAMED®

Ideal for vacuum mixing



PALAMED®	Description	Content	REF
PALAMED®	Medium-viscosity bone cement	1x40 2x40	66017786 66017784
PALAMED®G	Medium-viscosity bone cement with gentamicin	1x20 2x20 1x40 2x40	66017783 66017782 66017780 66017779

PALAMIX®	Description	Content	REF
PALAMIX® uno 1 x 40 PALAMIX® uno 1 x 80 PALAMIX® uno 1 x 120	Vacuum mixing system	10 10 10	66036742 66036743 66036744
PALAMIX® duo 2 x 40 PALAMIX® duo 2 x 80	Vacuum mixing system with 2 cartridges	10 10	66036745 66036746
PALAMIX® duo hip 1 x 40 / 1 x 80	Vacuum mixing system hip with 2 cartridges	10	66044411
Accessories	PALAMIX® cement gun PALAMIX® vacuum pump PALAMIX® nozzle medium PALAMIX® nozzle slim	1 1 10 10	66036163 66036748 66043960 66036747

Simply order from Heraeus.

www.heraeus-medical.com