

Safety Rails

Q. Can safety rails be used on ArjoHuntleigh systems?

A. ArjoHuntleigh recommends the use of safety rails on the bed when using their dynamic, active or static systems. Local facility policies regarding the use of safety rails should also be taken into account and their use considered in light of an individual risk assessment. Each facility should assess the fit of the device to the bed frame and side rails to ensure minimum risk of entrapment.

Note: that the side rails are used to indicate the edge of the bed and not to constrain or retain a confused or restless patient.



The list of FAQ's has been developed to answer common questions which arise during clinical practice and in the market place when clinical evidence is scarce. It is important to consider the holistic care of the patient and use clinical judgement making decisions based on the answers.

If you can not find an answer for which you are looking, please contact karen.milton@arjohuntleigh.com with your question(s) and we will endeavour to answer it/them and post it/them on this site to help educate others.