

ARJOHUNTLEIGH

GETINGE GROUP

CAREVO SHOWER TROLLEY



CAREVO™ shower trolley from ArjoHuntleigh is an innovative hygiene solution for immobile patients.

The new features, innovative design, and intuitive user interface provide an easy to use and efficient hygiene solution for both caregivers and immobile patients.

Designed with Emma (E) in mind.



...with people in mind

www.ArjoHuntleigh.com

KEY FEATURES

1 Integrated control panel



- Accessible at all times
- Simple, intuitive and efficient

2 Comfort handles



- Enhanced sense of security
- Promotes participation and mobility
- Provides support

3 Flexi zone



- Improved comfort
- Pressure relieving
- Easier handling and access

4 Ergo access area



- Work closer to your patient
- Reduce static load
- Improve work-flow

5 Elevated back support (Optional)



- Optimal patient comfort
- Easier and efficient patient transfers

6 Width adjustment option



- Greater versatility for more patient requirements
- Easier transfers and more handling options

7 Integrated shower handle holder



- Intuitive
- Easy to reach

8 Optimized head pillow and softer mattress



- Superior patient comfort
- Durable, easy to clean and disinfect

9 Simple daily maintenance



- Easy mattress attachment
- Easy mattress removal
- Easy to disinfect



Product Information

Max total load (SWL):	200 kg(440 lbs)
Max patient weight:	182 kg(401 lbs)
Weight (Carevo with mattress):	82 kg(181 lbs)
Width (stretcher):	586 mm(23")
Total width (incl. side supports):	759 mm(30")- 888 mm(35")
Total length:	2005 mm(79")
Length of chassis (incl. wheel frame):	1203 mm(47½")

Product Information

Width of chassis:	700 mm(27½")
Height of chassis (over wheels):	204 mm(8")
Height, top position floor to stretcher:	955 mm(37½")
Height, low position floor to stretcher:	605 mm(23¾")
Stroke range:	350 mm(13¾")
Drain hose holders included:	2 pieces
Max back support angle (Optional):	15°
Equipped with head pillow and foldable side supports with safety lock	

EXPLORE THE FULL RANGE OF ARJOHUNTLEIGH SHOWER SOLUTIONS AND SEE CAREVO IN ACTION

Visit www.ArjoHuntleigh.com/carevo



the preventable .com

Clinical research has shown that shower trolley solutions, like CAREVO, reduce the risk of serious injury for both caregivers and patients during hygiene tasks. For more information, visit the preventable.com

ARJOHUNTLEIGH

GETINGE GROUP

GETINGE GROUP is a leading global provider of products and systems that contribute to quality enhancement and cost efficiency within healthcare and life sciences. We operate under the three brands of ArjoHuntleigh, GETINGE and MAQUET. ArjoHuntleigh focuses on patient mobility and wound management solutions. GETINGE provides solutions for infection control within healthcare and contamination prevention within life sciences. MAQUET specializes in solutions, therapies and products for surgical interventions and intensive care.

Only ArjoHuntleigh designed parts, which are designed specifically for the purpose, should be used on the equipment and products supplied by ArjoHuntleigh. As our policy is one of continuous development we reserve the right to modify designs and specifications without prior notice.

® and ™ are trademarks belonging to the ArjoHuntleigh group of companies.
© ArjoHuntleigh, 2012.

www.ArjoHuntleigh.com



ArjoHuntleigh AB
Verkstadsvägen 5
SE-241 21 Eslöv
SWEDEN
Phone: +46 413 645 00
Fax: +46 413 645 63